ASLO Látex
 ASLO Latex / ASLO Látex
 Ref. 30.001.00 (KIT) / 30.001.00 (REAGENTE) **Responsável Técnico:**
 Dr. Gilson Serio Pizzo CRF MG – 5310
Anvisa 80027310268

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ENGLISH
BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE
 Kit for the determination of anti-streptolysin O (ASLO) antibodies in serum. *In vitro* diagnostic use.

STORAGE AND HANDLING
 • Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
 • Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
 • Do not use reagents whose shelf life has expired.
 • Ready-to-use reagents and controls.


WORKING PRINCIPLE
Method: Plate agglutination
 ASLO Latex is a plate agglutination method for qualitative and semi-quantitative detection, where the latex particles coated with streptolysin O are agglutinated by anti-streptolysin O antibodies present in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY
Sample Type: serum
Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

	Temperature	Stability Period
Serum	2 to 8 °C	7 days
	-20 °C	3 months

PRODUCT DESCRIPTION

 Suspension of Latex particles coated with Streptolysin O, pH 8.2. Preservative. The sensitivity of RLAX is traceable to NIBSC 97/662 (National Institute for Biological Standards and Control).
 Human serum with a concentration of ASO > 200 IU/ml. Preservative. (INCLUDED IN THE ASLO LATEX KIT)
 Human serum in buffer Phosphate, stabilizer and preservative. (INCLUDED IN THE ASLO LATEX KIT)

QUALITY CONTROL
 The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

CONTROL (+) (ASLO LÁTEX KIT)	 30.001.00
CONTROL (-) (ASLO LÁTEX KIT)	30.001.00

NECESSARY EQUIPMENT FOR TESTING
 • Reaction plate and disposable rod;
 • Pipette, clock or stopwatch..

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE
Automation: This procedure cannot be applied to automated analyzers.

Qualitative Method
 1. Before starting the tests, allow the RLAX, controls and samples to reach room temperature. Homogenize before use.
 2. On the reaction plate use a different circle for each sample or control.
 3. Place 20 µL of the sample / controls.
 4. In each circle add 20 µL of RLAX.
 5. Mix each determination using disposable rods (one for each sample).
 6. Homogenize the plate for 2 minutes in circular motions, tilting it back and forth. Observe the result under artificial light: Presence or absence of agglutination.

Semi-quantitative method

1. Dilute the sample with saline (0.9% NaCl) as follows:

Dilution	Saline (µL)	Sample(µL)	Dilution Factor
1/2	100	100 (pure)	2
1/4	100	100 (sol.1/2)	4
1/8	100	100 (sol.1/4)	8
1/16	100	100 (sol.1/8)	16

2) Carry out the test of each dilution according to the qualitative method.

B) CALCULATIONS
 The approximate concentration in the semi-quantitative test is obtained by multiplying the dilution factor by the sensitivity of the test (200 IU/mL).

Example:
 Aglutination up to 1/4 dilution.
 Dilution Factor = 4
 ASASLO = 800 IU / mL (4 x 200)

C) INTERPRETATION
 Streptolysin O is a toxic immunogenic exoenzyme produced by β-hemolytic Streptococcus groups A, C and G. The detection of ASLO antibodies is used for the diagnosis and treatment of diseases such as rheumatic fever, acute glomerulonephritis, and other streptococcal infections. These infections can later lead to impairment of the heart or kidneys. Early diagnosis, efficient treatment and patient monitoring can reduce these risks.

PERFORMANCE CHARACTERISTICS

Operating range
Up to 800 UI/mL

Sensitivity	Detection Limit
	200 UI/mL

Analytical Specificity			
Hemoglobin	Bilirubin	Triglycerides	Rheumatoid Factor
1000 mg/dL	40 mg/dL	2000 mg/dL	300 UI/mL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy

Number of Samples	100
Positive Agreement	98%
Negative Agreement	100%

Prozone High Dose Effect:
 There is no prozone effect up to 800 IU/mL of ASLO.

Precision:
Intra-assay: the study was carried out in two daily runs, in duplicate, for 10 days, using 04 samples, two positive and two negatives, with 100% agreement.
Inter-assay: the study was carried out in a daily run of four replicates, for 5 days, using 04 samples, two positive and two negatives, with 100% agreement

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS
 • Use protective equipment in accordance with the Good Laboratory Practices.
 • Follow the Good Laboratory Practices’ instructions to establish the quality of water.
 • Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
 • Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

	Adults	Up to 200 UI/mL
--	--------	-----------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS
 • Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
 • The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE
 • All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
 • To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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- BACH, G. L. et al. Latex Antistreptolysin O Test as a Tube Dilution Procedure. **Amer. J. Clin. Path.** v.57, p.209-211, 1972.
- HAFEEJEE. Rheumatic Fever and Rheumatic Heart Disease: The Current Status of its Immunology, Diagnostic Criteria, and Prophylaxis. **Quarterly Journal Of Medicine**, New Series 84, n° 305, p.643-658, 1992.
- KASSEM, Ahmed Samir; ZAHER, Salah Rafik. An International Comparison of the Prevalence of Streptococcal Infections and Rheumatic Fever in Children. **Pediatric Annals**. v. 21, n° 12, p.835-842, 1992.
- KLEIN, George C.; BAKER, Carolyn N.; JONES, Wallis L. "Upper Limits of Normal" Antistreptolysin O and Antideoxyribonuclease B Titers. **American Society for Microbiology**, v. 21, n° 6, p.999-1001, 1971.
- SPAUN, J. et al. International Standard for Antistreptolysin-O. **Bull. Wld Hlth Org.** v. 24, p.271-279, 1961.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: **AACC Press**, 2000.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES			
	Consultar as instruções para utilização Consult instructions for use Consúltense las instrucciones de uso		Descartar corretamente Dispose properly Desechar adecuadamente
REF	Número de catálogo Catalog number Número de catálogo		Reagente Reagent Reactivo
LOT	Código do lote Batch code Código de lote		Límite de temperatura Temperature limitation Límite de temperatura
	Risco biológico Biological risk Riesgo biológico		Nocivo / Irritante Harmful / Irritant Nocivo / Irritante
CONTROL	Controle Control Control		Validade Use by date Fecha de Caducidad
IVD	Produto para a saúde para diagnóstico <i>in vitro</i> In Vitro Diagnostic medical device Producto sanitario para diagnóstico <i>in vitro</i>		
	Atenção Attention Atención		