ASO Turbilátex

ASO Turbilátex / ASO Turbilátex
Ref. 20.004.00

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ENGLISH
 BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE
Kit for the determination of anti-streptolysin O (ASO) antibodies in serum. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- **Work Reagent:** mix in the proportion of 4 parts of RDIL + 1 part of RLAX. Homogenize gently. The work reagent is stable for 20 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).

WORKING PRINCIPLE
Method: Immunoturbidimetric
Latex particles coated with Streptolysin O are agglutinated by anti-streptolysin O (ASO) antibodies present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of ASO.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY
Sample Type: serum
Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

	Temperature	Stability Period
Serum	2 to 8 °C	7 days
	-20 °C	3 months

PRODUCT DESCRIPTION		
	Tris buffer 20 mmol / L pH 8.2, preservative.	
	Suspension of latex particles coated with Streptolysin O, pH 10,0, preservative.	
	Human serum containing anti-streptolysin O, preservative. The ASO value is traceable to the International Standard for anti-streptolysin O, NIBSC (97/662).	

QUALITY CONTROL
The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

ASO Turbilátex CAL (Kit Included)	20.004.00
Rheumatic Control Level 1	21.007.00
Rheumatic Control Level 2	21.008.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 540 nm (530 - 550 nm).

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE
Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica2da.ind.br

REACTION PARAMETERS	
Type of reaction	Fixed time
Wave length	540 nm
Temperature	37 °C
Sample volume/CAL*	3 µL
RT volume*	300 µL
Reading 1 (A1)	Immediately
Reading 2 (A2)	2 minutes

* Sample volumes/CAL and RT be modified as long as the established ratio is maintained.

B) CALCULATIONS
The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION
Streptolysin O is a toxic immunogenic exoenzyme produced by β-hemolytic *Streptococci* from groups A, C and G. The detection of ASLO antibodies is used for the diagnosis and

treatment of diseases such as rheumatic fever, acute glomerulonephritis, and other streptococcal infections. These infections can later lead to impairment of the heart or kidneys. Early diagnosis, efficient treatment, and patient monitoring can reduce these risks.

PERFORMANCE CHARACTERISTICS

Operating range	
34.2 to 736 UI/mL	

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
7,2 UI/mL	34,2 UI/mL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
1000 mg/dL	40 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 0.9999x - 0.053
Correlation Coefficient (R)	1.000

By applying the regression equation, the total systematic error estimated for decision levels of 120 UI/mL and 500 UI/mL was -0.15% and -0.12%, respectively.

Prozone High Dose Effect:
There is no prozone effect up to 4000 IU/mL of ASO.

Precision:
Determined with two runs in duplicate per day for 20 days.

Samples (UI/mL)	Replicates	Within-Run Precision		Total Precision	
		SD (UI/mL)	CV (%)	SD (UI/mL)	CV (%)
79.3	80	1.33	1.7	1.56	2.0
317	80	8.14	2.6	8.51	2.7
457	80	8.92	2.0	9.54	2.1

CV: Coefficient of variation; SD: Standard deviation

- RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**
- Use protective equipment in accordance with the Good Laboratory Practices.
 - Follow the Good Laboratory Practices' instructions to establish the quality of water.
 - Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
 - Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Adults	Up to 200 UI/mL
--------	-----------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

- WARNINGS AND PRECAUTIONS**
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
 - The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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- HAFEEJEE, I.; Rheumatic Fever and Rheumatic Heart Disease: The Current Status of its Immunology, Diagnostic Criteria, and Prophylaxis. *Quarterly Journal of Medicine.* v. 84, n. 305, p. 641-658, 1992.
- KLEIN, G.C.; *et al.* "Upper Limits of Normal" Antistreptolysin O and Antideoxyribonuclease B Titers. *Applied Microbiology.* v. 21, n. 6, p. 999-1001, 1971.
- SPAUN, J.; *et al.* International Standard for Antistreptolysin-O. *Bull Wld Hlth Org.* v. 24, p. 271-279, 1961.
- WANNAMAKER, L.W.; *et al.* Antibody Titers in Acute Rheumatic Fever. *Circulation.* v. 21, p. 598-614, 1960.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.

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