AST – TGO

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine the activity of Aspartate Transaminase (AST) or Glutamic-Oxalacetic Transaminase (GOT) in serum and plasma samples. *In vitro* diagnostic use.



STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 30 days.
- Do not use reagents whose shelf life has expired.
- Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 14 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

WORKING PRINCIPLE

Method: Kinetic UV

The Aspartate Transaminase or Glutamic-Oxaloacetic Transaminase (AST or GOT) catalyzes the transference of an amino group from aspartate to alpha-ketoglutarate to form oxaloacetate and glutamate. The enzymatic activity is determined by coupling this reaction to malate dehydrogenase (MDH) and is proportional to the consumption rate of NADH, which is spectrophotometrically measured at 340 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (EDTA and heparin)

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Serum and Plasma*	4 to 8 °C	7 days

*Samples collected with the citrate, sodium fluoride and sodium oxalate anticoagulants interfere in the dosage of AST.

PRODUCT DESCRIPTION

R 1	Tris buffer ≥ 80 mM, L-aspartic acid ≥ 100 mM; malate dehydrogenase ≥ 100 U/L; lactate dehydrogenase ≥ 500 U/L; alpha-ketoglutarate ≥ 10 mM; activators; stabilizers; surfactants; preservative.	 
	NADH ≥ 0,5 mM; stabilizer; preservative.	

The determination of AST-GOT activity is traceable to the reference material ERM-AD457/IFCC.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00	REF
Normal Control - Quantinorm	13.003.00	
Pathological Control - Quantialt	13.004.00	

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

- Heat the work reagent at 37 °C for 3 minutes.
- Pipette in the test tubes:

	Sample
Sample	100 µL
Work Reagent	1,0 mL

- Homogenize and insert it immediately in the thermostated cuvette at 37 °C.
- After 1 minute, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute ($\Delta A/\text{min}$):

$$\Delta A/\text{min} = \frac{(A0 - A1) + (A1 - A2) + (A2 - A3)}{3}$$

The GOT activity is calculate by multiplying the $\Delta A/\text{min}$ by the factor below:

$$\text{GOT (U/L)} = \Delta A/\text{min} \times 1746$$

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicauld.ind.br

C) INTERPRETATION

The Glutamate Oxaloacetate Transaminase (GOT) or Aspartate Transaminase (AST) is an intracellular enzyme present in high quantities in cardiac and skeletal muscles, liver, kidney and brain. Elevated levels of GOT help the diagnostic of cardiac, hepatic and muscle diseases. The activity of this enzyme in the myocardial infarction is elevated within the first 12 hours, reaching a peak in 24 hours and come back to normal in the fifth day. Discrete elevations are observed in pregnancy. Elevated levels are also found in hepatic necrosis, hemolytic anemias, acute pancreatitis, cirrhosis, hepatitis, obstructive jaundice, mononucleosis, hypothyroidism, brain necrosis and trauma, severe burns, muscle dystrophies, skeletal muscle injuries, cardiac caterisism and angioplasty. Several common drugs may increase the level of GOT (isoniazide, erythromycin, progesterone, anabolic steroids, etc.), being then used in monitoring therapies which uses hepatotoxic drugs.

PERFORMANCE CHARACTERISTICS

Operating range
3.05 to 440.0 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
2.41 U/L	3.05 U/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	40 mg/dL	350 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1.006x + 0.53$
Correlation Coefficient (R)	0.9996

By applying the regression equation, the total systematic error estimated for decision levels of 50 U/L and 150 U/L was 1.66% and 0.95%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
64.17	80	1.05	1.60	1.18	1.80
240.49	80	1.27	0.50	1.27	0.50
424.86	80	2.77	0.70	3.22	0.80

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	Normal Values (37 °C)
Men	< 35 U/L
Women	< 31 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

AST (U/L) x 0,017 = AST (µKat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- SCHUMANN, G., et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C. Part 5. Reference Procedure for the Measurement of Catalytic Concentration of Aspartate Amino-transferase. **Clin. Chem. Lab. Med.** v. 40 p.725-733, 2002.
- BURTSI, C. A.; ASHWOOD, E. R. Tietz textbook of clinical chemistry, 4 ed. Philadelphia: W.B.Saunders, 1994. 836p.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - vol. 2**, 5 ed. Washington DC: AAC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.

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