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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE


Kit intended to determine uric acid in serum and urine samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic – Trinder
Uricase catalyzes the oxidation of uric acid from sample, producing allantoin and hydrogen peroxide. The latter, in the presence of peroxidase, reacts with DHBS and 4-aminoantipyrine, producing a pink compound (quinonimine) that is spectrophotometrically measured at 505 nm. Color intensity is proportional to uric acid concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum and urine.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Serum	4 to 8 °C	5 days
	-20 °C	6 months
24 Hours Urine	20 a 25 °C	4 days

Preparation:
Urine: work with samples collected in the period of 24 hours. Homogenize the urine, measure its volume and separate 20 mL. Adjust the pH between 7 and 9 with 5% NaOH solution. In case of turbidity, incubate the sample for 10 minutes at 56 °C. Perform a 1:10 dilution of the aliquot with purified water and perform the assay. Multiply the result obtained by 10. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

PRODUCT DESCRIPTION

R 1	Phosphate buffer ≥ 50 mmol/L; uricase ≥ 100 U/L; peroxidase ≥ 1000 U/L; 4-aminoantipyrine ≥ 50 µmol/L; DHBS ≥ 0,5 mmol/L; activators; stabilizers and preservatives.	X
	Uric acid concentration of 10 mg/dL. Isopropilic alcohol ≥ 10% v/v; stabilizers and preservatives. Traceable to reference material NIST 913a.	X

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 505 nm (490 – 510 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	20 µL	-
Sample	-	-	20 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 10 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 505 nm (490-510 nm). The final reaction is stable for 15 minutes.

B) CALCULATIONS

Calculation for Serum

Uric Acid (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:
Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Uric Acid (mg/dL) = Sample Absorbance x Concentration Factor

Calculations for urine samples:

Uric Acid (mg/24 hours) = $\frac{\text{Uric Acid (mg/dL)} \times \text{Urinary Volume (mL)}}{100}$

* Determine the concentration of uric acid in the urine sample in mg/dL using the same calculation procedure for serum. Multiply the result by 10.

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnica2tda.ind.br

C) INTERPRETATION

Uric acid is the product of purine metabolism in human organism. Its dosage is used predominantly to investigate gout, a clinical condition that results from a primary hyperuricemia or is caused by other conditions or treatments which result in secondary hyperuricemia. Furthermore, uric acid determination is used in the diagnosis and monitoring of numerous renal and metabolic disorders, such as: renal failure, leukemia, psoriasis, malnutrition, in patients using cytotoxic drugs and in pregnancy-induced hypertension (pre-eclampsia).

PERFORMANCE CHARACTERISTICS

	Operating range
Serum	0.23 to 20.00 mg/dL
Urine	9.19 to 210.10 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

	Sensitivity	
	Detection Limit	Quantification Limit
Serum	0.18 mg/dL	0.23 mg/dL
Urine	0.146 mg/dL	9.19 mg/dL

	Analytical Specificity		
	Hemoglobin	Bilirubin	Triglycerides
Serum	200 mg/dL	6 mg/dL	900 mg/dL
Urine	200 mg/dL	15 mg/dL	-

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine
Number of Samples	40 in duplicate	27 in duplicate
Regression Equation	$y = 0.991x + 0.007$	$y = 0.987x - 1.762$
Correlation Coefficient (R)	0.9901	0.9944

Serum: using the regression equation obtained, the total systematic errors estimated for decision levels of 8.0 mg/dL and 15.0 mg/dL were, respectively, -0.81% and -0.85%.
Urine: using the regression equation obtained, the total systematic errors for decision levels of 300.0 mg/24h and 650.0 mg/24h were, respectively, -1.89% and -1.57%.

Precision:

The study was carried out in two runs per day, in duplicate, for 20 days.

Serum Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
3.89	80	0.026	0.70	0.057	1.50
6.39	80	0.054	0.80	0.090	1.40
12.51	80	0.070	0.60	0.185	1.50

Urine Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
13.26	80	0.265	2.0	0.332	2.2
19.41	80	0.498	2.6	0.515	2.7

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	Serum	Urine
	Male	3.5 – 7.2 mg/dL 208 – 428 mmol/L
Female	2.6 – 6.0 mg/dL 155 – 357 mmol/L	

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Uric Acid (mg/dL) X 59.5 = Uric Acid (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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