Ácidos Biliares Totais (TBA) VET

Total Bile Acids (TBA) VET/Ácidos Biliares Totales (TBA) VET
Ref. 90.055.00

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Caninos	< 5,0 µmol/L
Felinos	< 5,0 µmol/L
Equinos	5,0 -28,0 µmol/L
Bovinos	20,0 – 80,0 µmol/L

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ENGLISH

INTENDED USE

Kit intended to determine total bile acids (TBA) in serum and plasma samples. Veterinary diagnostic use only


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable for 4 weeks, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic

The enzyme 3-α-hydroxysterol dehydrogenase (3α-HSD), in the presence of Thio-NAD, oxidize the bile acids from sample to form 3-ketosterol and Thio-NADH. The Thio-NADH rate of formation, which is spectrophotometrically measured at 405 nm, is proportional to the bile acids concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma EDTA.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material. Samples should be collected while fasting (8 hours), as the bile acids concentration increases after meals.
Preservation:

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	7 days
	-20 °C	1 month

PRODUCT DESCRIPTION

R 1	Goods buffer ≥ 5 mmol/L; Thio-NAD ≥ 900 mg/dL; stabilizers; activators; preservatives.	X
R 2	Goods buffer ≥ 300 mmol/L; Thio-NADH ≥ 6 mg/dL; 3α-HSD ≥ 12.000 U/L; stabilizers; activators; preservatives.	X
STD	Bile acids; stabilizers; preservatives. The STD concentration is printed on the label.	X

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. It is recommended to perform the internal laboratory quality control to ensure the obtained results.

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 and 660 nm.
- Glass and/or automatic pipettes.
- Test tubes.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaltda.ind.br

Reaction type	Fixed Time / Positive
Wavelength	405 / 600 nm
Temperature	37 °C
Sample/STD volume*	3 µL
R1 volume*	225 µL
Sample/STD + R1 Incubation Time	300 seconds
R2 Volume*	75 µL
First Reading	300 seconds
Second Reading	120 seconds
Calibration	One Point

* Volumes can be modified as long as the established proportion is maintained.

B) CALCULATIONS

Total Bile Acids (µmol/L) = $\frac{\text{Sample}'(A2 - A1)}{\text{STD}'(A2 - A1)}$ x Concentration of STD (µmol/L)

Calculations with de Calibration Factor:

Calibration Factor = $\frac{\text{Concentration of STD (µmol/L)}}{\text{STD}'(A2 - A1)}$

Total Bile Acids (µmol/L) = *Calibration Factor* x *STD'*(A2 – A1)

C) INTERPRETATION

Total Bile Acids (TBA) are chemical compounds derived from cholesterol and synthesized in the liver. They play a fundamental role in animal metabolism, especially in the process of digesting and absorbing fats. TBAs are secreted in bile and act as emulsifying agents, aiding in the breakdown and solubilization of dietary fats. Additionally, they also play a role in the absorption of fat-soluble vitamins such as vitamins A, D, E, and K. Total Bile Acids (TBA) are sensitive indicators of liver function and biliary system health. This means they can be used to detect early disorders like acute hepatitis, chronic hepatitis, liver cirrhosis, steatosis, liver cancer, and biliary obstructions. Furthermore, monitoring TBAs is particularly valuable in clinical situations where there is suspicion of liver impairment, such as in cases of poisoning, metabolic diseases, and certain infections. Despite total bile acids being considered sensitive and specific markers for liver disorders, the information obtained from this test should be confirmed by other liver function tests for diagnosis confirmation.

PERFORMANCE CHARACTERISTICS

Operating range
1.58 a 194.78 µmol/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
1.50 µmol/L

Analytical Specificity

Species	Hemoglobin	Bilirubin	Triglycerides
Canines	500 mg/dL	10 mg/dL	750 mg/dL
Felines	500 mg/dL	10 mg/dL	1000 mg/dL
Equines	500 mg/dL	10 mg/dL	750 mg/dL
Bovines	500 mg/dL	40 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy		
Control	External Control 1	External Control 2
Interval (µmol/L)	20.2 – 30.2 (média 25.2)	33.1 – 49.7 (média 41.4)
Average (µmol/L)	24.05	40.89
%DIF	-4.6%	-1.2%

%DIF less than the total error allowed for the analyte: 10%.

Accuracy		
Number of Samples	40 in duplicate	
Regression Equation	y = 1.021x + 0.135	
Correlation Coefficient (R)	0.999	

Precision:

The studies were performed in two runs, in duplicate.

Samples (µmol/L)	Within-Run Precision		Samples (µmol/L)	Between-Run Precision	
	SD	CV (%)		SD	CV (%)
4.91	0.15	3.0	4.86	0.2	4.1
70.81	1.26	1.9	71.57	1.98	2.8

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.

- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.
- Reagent contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
- Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide from building up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.

REFERENCE RANGES

Canines	< 5,0 µmol/L
Felines	< 5,0 µmol/L
Equines	5,0 -28,0 µmol/L
Bovines	20,0 – 80,0 µmol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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- BURTIS, C. A.; ASHWOOD, E. R.; BRUNS, D. E. **Tietz Fundamentos de Química Clínica**, Saunders Elsevier, 6 ed, 2008.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACCPress, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.
- LIU, Yanan et al. Detection technologies and metabolic profiling of bile acids: a comprehensive review. Lipids in health and disease, v. 17, n. 1, p. 1-13, 2018.
- DANESE, Elisa et al. Bile acids quantification by liquid chromatography–tandem mass spectrometry: Method validation, reference range, and interference study. Diagnostics, v. 10, n. 7, p. 462, 2020.
- ZHANG, Guo-Hua et al. An enzymatic cycling method for the determination of serum total bile acids with recombinant 3α-hydroxysteroid dehydrogenase. Biochemical and biophysical research communications, v. 326, n. 1, p. 87-92, 2004.