## Albumina VET

Albumin VET / Albumina VET  
Ref. 90.020.00

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#### ENGLISH

##### INTENDED USE

Kit intended to determine albumin in serum samples. Veterinary diagnostic use only.

##### STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30 °C) are followed.
- Do not use reagents whose shelf life has expired.

##### WORKING PRINCIPLE

**Method:** Bromocresol Green

In acidic pH, albumin from sample reacts with bromocresol green to form a blue-green complex that is spectrophotometrically measured at 630 nm. The intensity of the color is proportional to the albumin concentration in the sample.

##### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum.

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum	4 to 8 °C	3 days
	-20 °C	7 days

#### PRODUCT DESCRIPTION

**R 1**

Succinate buffer ≥ 20 mmol/L; bromocresol green ≥ 100 µmol/L; stabilizer; detergent; preservative.

**STD**

Phosphate buffer ≥ 20 mmol/L; bovine albumin at a concentration equivalent to 4,0 g/dL; preservative. Traceable to reference material ERM-DA470k/IFCC.

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal Vet	<b>REF</b>	90.039.00
Normal Control – Quantinorm Vet		90.040.00
Pathological Control – Quantialt Vet		90.041.00

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 630 nm (600 – 640 nm).
- Test tubes, glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	5 µL	-
Sample	-	-	5 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 15 to 30 °C for 10 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 630 nm (600 – 640 nm). The final reaction is stable for 30 minutes.

#### B) CALCULATIONS

Albumin (g/dL) =  $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (g/dL)}$

**Calculations with the Calibration Factor:**

Calibration Factor =  $\frac{\text{STD Concentration (g/dL)}}{\text{STD Absorbance}}$

Albumin (g/dL) = Sample Absorbance x Calibration Factor

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

#### C) INTERPRETATION

It is the most abundant plasma protein, synthesized mainly in the liver. Its main function is to maintain colloidal osmotic pressure inside blood vessels and extravascular spaces. Albumin transports many compounds, such as: free fatty acids, phospholipids, metal ions, amino acids, drugs, hormones and bilirubin. Elevated albumin concentrations occur only in acute dehydration. Its reduction occurs in several clinical conditions such as: analbuminemia, inflammation, liver disease, urinary and gastrointestinal losses, malnutrition, edema, ascites. Its low concentration in the blood has an important pharmacokinetic effect, since Albumin can bind to numerous drugs.

#### PERFORMANCE CHARACTERISTICS

Operating range
0.22 to 4.98 g/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
0.22 g/dL

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	500 mg/dL	40 mg/dL	1250 mg/dL
Felines	500 mg/dL	40 mg/dL	2000 mg/dL
Equines	500 mg/dL	40 mg/dL	750 mg/dL
Bovines	300 mg/dL	40 mg/dL	750 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	40 in duplicate
Regression Equation	$y = 0.9521x + 0.1858$
Correlation Coefficient (R)	0.9750
Accuracy – Felines	
Number of Samples	20 in duplicate
Regression Equation	$y = 0.9943x + 0.0156$
Correlation Coefficient (R)	0.9975
Accuracy – Equines	
Number of Samples	11 in duplicate
Regression Equation	$y = 1.035x - 0.1152$
Correlation Coefficient (R)	0.9938
Accuracy - Bovines	
Number of Samples	10 in duplicate
Regression Equation	$y = 0.9833x + 0.0411$
Correlation Coefficient (R)	0.9976

**Precision:**

Determined with two runs in 20 replicates.

Samples (U/L)	Number of samples	Within-Run Precision		Samples (U/L)	Number of samples	Total Precision	
		SD (U/L)	CV (%)			SD (U/L)	CV (%)
2,56	20	0,05	1,80	2,50	40	0,08	3,10

4.76	20	0.13	2.70	4.73	40	0.16	3.40
------	----	------	------	------	----	------	------

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

#### REFERENCE RANGES

Canine	2.6 – 3.3 g/dL
Feline	2.1 – 3.3 g/dL
Bovine	3.03 – 3.55 g/dL
Equine	2.6 – 3.70 g/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units (g/L):

Albumin (g/dL) X 10 = Albumin (g/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnicaltda.com.br](mailto:sac@biotecnicaltda.com.br).

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