Albumina

Albumin / Albúmina
Ref. 10.002.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE
Kit intended to determine albumin in serum samples. Diagnostic use only.
STORAGE AND HANDLING
Store at 15 to 30 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
Reagent ready for use.
Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30 °C) are followed.
Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE
Method: Bromocresol Green
Albumin reacts with bromocresol green at acidic pH to form a blue-green complex that is spectrophotometrically measured at 630 nm. The intensity of the color is proportional to the concentration of albumin in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY
Sample Type: serum.
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

Table: Stability Period (Serum: 4 to 8 °C: 3 days; -20 °C: 7 days)

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum.
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

Table: Stability Period (Serum: 4 to 8 °C: 3 days; -20 °C: 7 days)

PRODUCT DESCRIPTION

R 1: Succinate buffer ≥ 20 mmol/L; bromocresol green ≥ 100 µmol/L; stabilizer; detergent; preservative.
STD: Phosphate Buffer ≥ 20 mmol/L; bovine albumin at a concentration equivalent to 4,0 g/dL; preservative. Traceable to reference material ERM-DA470k / IFCC.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Table: Control Name, REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 630 nm (600 – 640 nm).
Test tubes.
Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

- 1. Pipette in the test tubes:
Table with 4 columns: STD, Sample, Blank, Standard
2. Homogenize and incubate at 15 to 30 °C for 10 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 630 nm (600-640 nm). The final reaction is stable for 30 minutes.

B) CALCULATIONS

Albumin (g/dL) = Sample Absorbance x STD Concentration (g/dL) / STD Absorbance

Calculations with the Calibration Factor:

Calibration Factor = STD Concentration (g/dL) / STD Absorbance

Albumin (g/dL) = Sample Absorbance x Concentration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechincaldta.ind.br

C) INTERPRETATION

Albumin is the most abundant plasma protein, synthesized primarily in the liver. Its main function is the maintenance of the coloidal osmotic pressure inside the blood vessels and extravascular spaces. It also carries several compounds, such as: free fatty acids, phospholipids, metal ions, amino acids, drugs, hormones and bilirubin. Elevated albumin concentrations occur only in acute dehydration and is of no clinical importance. Its reduction occurs in several clinical conditions, such as: analbuminemia, inflammation, liver disease, urinary and gastrointestinal losses, malnutrition, edema and ascites. Albumin can bind to numerous drugs, which is why low blood albumin concentrations have a significant pharmacokinetic effect.

PERFORMANCE CHARACTERISTICS

Table: Operating range (0.27 to 6.00 g/dL)

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Table: Sensitivity (Detection Limit: 0.06 g/dL; Quantification Limit: 0.27 g/dL)

Table: Analytical Specificity (Hemoglobin: 500 mg/dL; Bilirubin: 40 mg/dL; Triglycerides: 2000 mg/dL)

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Table: Accuracy (Number of Samples: 40 in duplicate; Regression Equation: y = 1.015x - 0.057; Correlation Coefficient (R): 0.9979)

By applying the regression equation, the total systematic error estimated for decision levels of 2,0 g/dL and 4,0 g/dL was -1.35% and 0.08%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Table: Precision data (Samples, Replicates, Within-Run Precision, Total Precision)

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
Follow the Good Laboratory Practices' instructions to establish the quality of water.
Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Table: Reference Range (Serum: 3.5 - 5.2 g/dL)

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (g/L):

Albumin (g/dL) X 10 = Albumin (g/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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