## Alfa Amilase VET

Alpha Amylase VET / Alfa Amilase VET  
 Ref. 90.018.00

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$$\text{Alfa-Amilase (U/L)} \times 0,017 = \text{Alfa-Amilase (µKat/L)}$$

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal Vet	90.039.00
Normal Control - Quantinorm Vet	90.040.00
Pathological Control - Quantilat Vet	90.041.00

REF

90.039.00  
90.040.00  
90.041.00

#### NECESSARY EQUIPMENT FOR TESTING

##### A) TEST PROCEDURE

1. Heat the reagent (R1) at 37 °C for 3 minutes.
2. Pipette in the test tubes:

Sample
10 µL
1.0 mL

3. Homogenize and insert into the cuvette holder thermostated at 37 °C, zeroed with purified water at 405nm. Set the stopwatch.

4. After 1 minute, measure the initial absorbance A<sub>0</sub> and read it again every minute, for 3 minutes (A<sub>1</sub>, A<sub>2</sub> and A<sub>3</sub>, respectively).

#### B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute ( $\Delta A/min$ ):

$$\Delta A/min = (A_1 - A_0) + (A_2 - A_1) + (A_3 - A_2)$$

3

The Alpha-Amylase activity is calculated by multiplying the  $\Delta A/min$  by the factor below:

$$\text{Alpha-Amilase (U/L)} = \Delta A/min \times 6520$$

#### C) INTERPRETATION

Alpha-amylase is an enzyme of the hydrolases class that catalyzes the hydrolysis of 1,4- $\alpha$ -glycosidic bonds in polysaccharides. The amylase activity present in serum and urine has pancreatic and salivary gland origins, increasing greatly in acute pancreatitis and inflammation of this gland. Increased values are also found in several acute intra-abdominal disorders, such as intestinal obstruction, renal failure, urinary obstruction, and in several extrahepatic conditions and in neoplastic diseases. Some drugs can cause pancreatitis, and consequent hyperamylasemia. Amylase alteration is not a good means of evaluating pancreatitis in cats. In this species, complementary exams should be used.

#### PERFORMANCE CHARACTERISTICS

Operating range
3.87 to 1513.80 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit

$$3.85 \text{ U/L}$$

Analytical Specificity
Species
Canines
Felines
Equines
Bovines

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines

Accuracy - Equines	
Number of Samples	11 em duplicata
Regression Equation	y = 1.0057x + 1.5388
Correlation Coefficient (R)	0.9991
Accuracy - Bovines	
Number of Samples	10 em duplicata
Regression Equation	y = 0.9693x + 2.1735
Correlation Coefficient (R)	0.9970

#### Precision:

Determined with two runs in 20 replicates.

Samples (U/L)	Number of samples	Within-Run Precision		Samples (U/L)	Number of samples	Total Precision	
		SD (U/L)	CV (%)			SD (U/L)	CV (%)
236,87	20	3,11	1,30	235,41	40	3,73	1,60
1102,02	20	41,10	3,70	1090,67	40	36,99	3,40

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Canine	185 - 700 U/L
Feline	< 500 U/L
Bovine	< 31 U/L
Equine	75 - 150 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units ( $\mu$ Kat/L):

Alpha-amylase (U/L)  $\times 0,017$  = Alpha-amylase ( $\mu$ Kat/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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