Alfa Amilase

Alpha Amylase / Alfa Amilosa
 Ref. 11.001.00 / FD38.047.00 / SD39.056.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine the activity of Alpha-Amylase in serum, plasma and urine samples. *In vitro* diagnostic use.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable for 30 days, as long as the recommended storage conditions (2 to 8 °C) are followed.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 28 days.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Kinetic

Alpha-amylase catalyzes the hydrolysis of α-(2-chloro-4-nitrophenyl)-β-1,4-galactopyranosylmaltoside (GalG2-CNP) into 2-chloro-4-nitrophenyl (CNP) and 1,4-galactopyranosylmaltoside (GalG2). The enzymatic activity is determined by the formation rate of CNP, which is spectrophotometrically measured at 405 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum, heparin plasma and urine

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Serum and Plasma	4 to 8 °C	7 days
	20 to 25 °C	7 days

Urine	-20 °C	1 year
	4 a 8 °C	10 days
	20 to 25 °C	2 days
	-20 °C	3 weeks

PRODUCT DESCRIPTION

R 1	Buffer MES > 20 mM; potassium thiocyanate > 100mM; GalG2-CNP > 2 mM; activators; surfactants; preservative	X

The determination of Alpha-amylase activity is traceable to the specific absorptivity of the CNP chromogen using calibrated pipettes and manual spectrophotometer that provide absolute values.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Serum and Plasma	Urine
R1	1,0 mL	1,0 mL
Sample	20 µL	10 µL

- Homogenize and insert into the cuvette holder thermostated at 37 °C, zeroed with purified water at 405nm. Set the stopwatch.
- After 1 minute, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

The Alpha-Amylase activity is calculated by multiplying the ΔA/min by the factor below:

	Amylase (U/L)
Serum and Plasma	ΔA/min x 3292
Urine	ΔA/min x 6520

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Alpha-amylase is a hydrolases-class enzyme that catalyzes the hydrolysis of 1,4-α-glycosidic bonds in polysaccharides. Present in serum and urine has both pancreatic and salivary gland origins. Increased values are found in acute pancreatitis and inflammation of this gland, beyond several acute intra-abdominal disorders, in various extrahepatic conditions and in neoplastic diseases.

PERFORMANCE CHARACTERISTICS

	Operating range
Serum and Plasma	4,20 a 1200 U/L
Urine	5,20 a 2000 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

	Sensitivity	
	Detection Limit	Quantification Limit
Serum and Plasma	2,10 U/L	4,20 U/L
Urine	1,26 U/L	5,20 U/L

	Analytical Specificity		
	Hemoglobin	Bilirubin	Triglycerides
Serum and Plasma	200 mg/dL	40 mg/dL	1250 mg/dL
Urine	500 mg/dL	40 mg/dL	-

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum and Plasma	Urine
Number of Samples	40 in duplicate	40 in duplicate
Regression Equation	y = 0.992x - 1.702	y = 0.989x - 4.316
Correlation Coefficient (R)	0.9916	0.9980

Serum and Plasma: by applying the regression equation, the total systematic error estimated for decision levels of 40 U/L and 100 U/L was -5.05% and -2.50%, respectively. **Urine:** by applying the regression equation, the total systematic error estimated for decision levels of 100 U/L and 300 U/L was -5.42% and -2.54%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Serum and Plasma Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
114.60	80	0.555	0.50	1.303	1.10
214.90	80	1.209	0.60	2.690	1.30
495.20	80	0.971	0.20	2.987	0.60

Urine Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
217.70	80	2.181	1.00	2.856	1.30
430.60	80	4.133	1.00	5.507	1.30
832.70	80	7.233	0.90	8.628	1.00

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- From different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated, and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.
- Additional care should be taken to avoid R1 contact with saliva or sweat.
- Do not mix reagents from different lots or exchange the caps specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Serum and Plasma	Urine
< 100 U/L	< 491 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

Alpha-amylase (U/L) x 0.017 = Alpha-amylase (µKat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and all others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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