Alfa Hidroxibutirato Desidrogenase

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine the activity of Alpha Hydroxybutyrate Dehydrogenase in serum samples. Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic

Alpha Hydroxybutyrate dehydrogenase (α-HBDH) catalyzes the reduction of α-ketobutyrate to α-hydroxybutyrate with concomitant oxidation of NADH to NAD⁺. The consumption rate of NADH, which is spectrophotometrically measured at 340 nm is proportional to the activity of α-HBDH in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	2 to 8 °C	7 days
	-20 °C	42 days

PRODUCT DESCRIPTION

R 1	α-ketobutyrate 5,6 mmol/L; phosphate 57 mmol/L; preservative 1,0 g/L.
R 2	NADPH I 10 mmol/L.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. It is recommended to perform the internal laboratory quality control to ensure the obtained results.

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnica.ind.br

Reaction Type	Kinetic / Negative
Wavelength	340 nm
Temperature	37 °C
Sample volume*	10 µL
R1 Volume*	200 µL
Sample + R1 Incubation Time	210 seconds
R2 Volume*	40 µL
R1 + Sample + R2 Incubation Time	30 seconds
Measure the absorbance continuously for 120 seconds	
Calibration	2 points

* Volumes can be modified as long as the established proportion is maintained.

Factor calculation:

For calibration, use the k factor (theoretical k = 4019). In different test conditions the k value can be changed. It is recommended to calculate the factor.

$$Factor = \frac{Total\ volume\ (mL) \times 10^7}{Sample\ volume\ (mL) \times \epsilon \times C\ (cm)}$$

Where:

Total volume = R1 + R2 + sample

ϵ = NADH absorptivity

C = optical path length of the cuvette

Exemplo:

$$Fator = \frac{0,250\ (mL) \times 10^7}{0,10\ (mL) \times 6220 \times 0,5} = 8038,58$$

B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

C) INTERPRETATION

Alpha Hydroxybutyrate Dehydrogenase (α-HBDH) is a Lactate Dehydrogenase (LDH) isoenzyme that has a higher affinity for α-ketobutyrate as substrate. Comparing with other LDH isoenzymes, α-HBDH is found at higher concentration in the cardiac muscle, while the other isoenzymes are found in the liver. Therefore, α-HBDH is a more specific marker for the diagnostic of myocardial infarction. Cardiac and hepatic diseases can be differentiated by calculating the α-HBDH/LDH ratio, in which decreased results indicate hepatic disease and increases results indicate cardiac disease. α-HBDH alongside LDH, Aspartate Transaminase (AST), Creatine Kinase (CK) and Creatine Kinase B Subunit (CK-MB) constitute the myocardial enzyme spectrum.

PERFORMANCE CHARACTERISTICS

Operating range
20 to 750 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Detection Limit
72 U/L

Analytical Specificity		
Bilirubin	Hemoglobin	Triglycerides
37,5 mg/dL	Any concentration	300 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y = 0.9596x + 4.6214
Correlation Coefficient (R)	0.9932

Precision:

Within-run precision was determined with 3 lots and 10 replicates in an analytical run; between-run precision was determined with 3 lots and 3 replicates in an analytical run.

Samples (U/L)	Within-Run Precision		Between-Run Precision	
	SD (U/L)	CV (%)	SD (U/L)	CV (%)
410.1	1.4	0.35	1.17	0.58

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Soro	72 – 182 U/L
------	--------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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