Autocal H

Autocal H / Autocal H

Ref. 13.002.00

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ENGLISH

INTENDED USE

Kit intended for calibration of quantitative biochemical assays. Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- After reconstituted, the analytes are stable for:
Bilirubins and Alkaline Phosphatase: 48 hours when stored at 2 to 8 °C, in the original flask.
Other analytes: 5 days when stored at 2 to 8 °C, in the original flask, and 15 days when stores at -20 °C, in microtubes.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Human lyophilized serum used to calibrate biochemical assays.

PREPARATION AND HANDLING OF THE PRODUCT

Reconstitute with 5 mL of purified water and let stand for 30 minutes. Gently mix by inversion, avoiding foam formation, until complete dissolution of the calibrator is observed.

PRODUCT DESCRIPTION

CAL

Human serum plus tissue extracts of human/animal origin.

The calibrator shows negative results for the HBs antigen and for the HIV and HCV antibodies. However, they should be treated as potentially infectious. Calibrator values are traceable to reference methods and reference materials from the National Institute of Standard and Technology (NIST) and European Reference Materials (ERM).

NECESSARY EQUIPMENT FOR TESTING

- Glass and/or automatic pipettes.
- Purified water.

TEST PROCEDURE

The calibrator should be used in accordance with the instructions for use of the corresponding biochemical method.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the controls below:

Controle Normal – QUANTINORM	REF	13.003.00
Controle Patológico – QUANTIALT		13.004.00

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix calibrators from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the product if it displays any signs in disagreement with the ones specified in its MSDS.
- Avoid leaving calibrators outside the specified storage conditions.

WARNINGS AND PRECAUTIONS

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646.
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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- CEMBROWSKI, G. Thoughts on quality-control systems: A laboratorians perspective. Clin. Chem. v.43, p.886-892, 1997.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981. Approved Guideline. CLSI document H47-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2008

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