Autocal VET

Autocal VET / Autocal VET
Ref. 90.039.00

FINALIDADE

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ENGLISH

INTENDED USE

Kit intended for calibration of quantitative biochemical assays. Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- After reconstituted, the analytes are stable for: Alkaline Phosphatase: 48 hours when stored at 2 to 8°C, in the original flask. Other analytes: 5 days when stored at 2 to 8°C, in the original flask, and 15 days when stores at -20 °C, in microtubes. Once defrosted, the analytes are stable for 24 hours.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Human lyophilized serum used to calibrate biochemical assays.

PREPARATION AND HANDLING OF THE PRODUCT

Leave the calibrator in temperature of 10 to 30 °C for 10 minutes. Reconstitute with 5.0 mL of purified water. Homogenize gently by inversion, avoiding foaming. Wait 30 minutes until complete dissolution of the product.

PRODUCT DESCRIPTION

CAL

Human serum plus tissue extracts of human/animal origin.

The calibrator shows negative results for the HBs antigen and for the HIV and HCV antibodies. However, they should be treated as potentially infectious. Calibrator values are traceable to reference methods and reference materials from the National Institute of Standard and Technology (NIST) and European Reference Materials (ERM).

NECESSARY EQUIPMENT FOR TESTING

- Glass and/or automatic pipettes.
- Purified water.

TEST PROCEDURE

The calibrator should be used in the same way as patient samples, according to the procedure established for the assays and in accordance with the equipment protocol.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the controls below:

Normal Control – Quantinorm Vet	REF	90.040.00
Pathological Control – Quantialt Vet		90.041.00

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix controls from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the product if it displays any signs in disagreement with the ones specified in its MSDS.
- Avoid leaving controls outside the specified storage conditions.

WARNINGS AND PRECAUTIONS

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646.
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper lace for potentially infectious material.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicaltda.com.br.

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