Bicarbonato - CO₂

Bicarbonato-CO₂ / Bicarbonato-CO₂
Ref. 12.016.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine Bicarbonate (CO₂) in serum and plasma samples. Diagnostic use only.



STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- It is recommended to calibrate the product daily.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic – Phosphoenolpyruvate carboxylase (PEPC)

Bicarbonato, in the presence of phosphoenolpyruvate and under catalytic action of phosphoenolpyruvate carboxylase (PEPC), is converted to oxaloacetate. Malate dehydrogenase (MDH) catalyzes the reduction of oxaloacetate to malate, oxidizing NADH to NAD⁺ in the process. The bicarbonate concentration in the sample is determined by the consumption rate of NADH, which is spectrophotometrically determined at 405 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma heparin

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	2 to 8 °C	1 week
	-20 °C	2 weeks

It is recommended to perform the Bicarbonate assay within 2 hours after the sample is collected. Avoid leaving the collection tube open, as exposure to air can cause pre-analytical problems.

PRODUCT DESCRIPTION

R 1	Tris-HCl buffer 25 mmol/L pH 7,5; Phosphoenolpyruvate 12,5 mmol/L; NADH 0,6 mmol/L; Malate dehydrogenase < 4100 U/L; Phosphoenolpyruvate carboxylase 400 U/L
------------	--

CAL	Hepes buffer 50 mmol/L; NaHCO ₃ 20-30 mmol/L. Traceable to reference material GBW(E)060023 do CNRM.
------------	--

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the control below:

Calibrator - CO ₂ CAL	REF	12.016.00
----------------------------------	------------	-----------

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnica.com.br

Reaction Type	Kinetic / Negative
Wavelength	405 nm
Temperature	37 °C
Sample volume*	2 µL
R1 Volume*	200 µL
Add R1 + Sample, measure the initial absorbance (A1) and after 60 seconds (A2)	
Calibration	2 points

* Volumes can be modified as long as the established proportion is maintained.

B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

C) INTERPRETATION

Bicarbonato (HCO₃⁻) is an electrolyte that takes part in the main buffering system in the blood, the carbonic acid-bicarbonate system. Bicarbonate's determination reflects the total concentration of CO₂ in the blood, so that it is used to evaluate the acid-base balance (pH) in the organism. Increased concentrations of CO₂ in blood (hypercapnia) might be due to pulmonary diseases, CO₂ inhalation, respiratory acidosis or metabolic alkalosis. Alongside with sodium, potassium and chloride, bicarbonate maintains the electric neutrality at a cellular level. Its determination aids in the diagnostics of electrolytic imbalance, acidosis or alkalosis. Additional assays allow to differentiate between respiratory or metabolic diseases.

PERFORMANCE CHARACTERISTICS

Operating range
5,3 a 50,0 mmol/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1,0 mmol/L	2,74 mmol/L

Analytical Specificity			
Ascorbic Acid	Hemoglobin	Triglycerides	Bilirubin
50 mg/dL	600 mg/dL	1000 mg/dL	50 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y=1,0102x - 0,4491
Correlation Coefficient (R)	0,9981

The product's accuracy was evaluated with the reference material GBW(E)060023 CNRM, at the decision level of 21.0 mmol/L, with a total systematic error of -2.70%.

Precision:

Within-run precision was determined with 20 replicates in an analytical run; between-run precision was determined with 5 replicates for 5 days.

Samples (mmol/L)	Within-Run Precision		Between-Run Precision	
	SD (mmol/L)	CV (%)	SD (mmol/L)	CV (%)
14.30	1.02	7.30	0.637	4.50
20.50	1.18	5.80	0.787	4.00

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Serum and Plasma	22 - 29 mmol/L
------------------	----------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the Unit mEq/L:

Bicarbonato (mmol/L) = Bicarbonato (mEq/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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