Bilirrubina Total

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ENGLISH

⚠️ BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine total bilirubin in serum and plasma (heparin and EDTA) samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: DCA (Diazotized Dichloroaniline).

In acidic pH, direct and indirect bilirubin are coupled with diazotized dichloroaniline, being solubilized and producing a colored complex, called azobilirubin, that can be measured at 550 nm (505-580 nm). The color intensity is proportional to the concentration of total bilirubin in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (heparin and EDTA).

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	7 days
	-20 °C	6 months

*Necessary to keep out of the light when stored for longer than 8 hours.

PRODUCT DESCRIPTION

R 1	Surfactants, stabilizers, catalyst and detergent
------------	--

R 2	Diazotized dichloroaniline < 2,0 mmol/L and stabilizer. pH < 5
------------	--

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal H	REF	13.002.00
Normal Control - Qualitnorm		13.003.00
Pathological Control - Quantitalt		13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 550 nm (505 – 580 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Sample Blank	Sample	Calibrator Blank	Calibrator
Calibrador	---	---	---	50 µL	50 µL
Sample	---	50 µL	50 µL	---	---
Purified water	50 µL	---	---	---	---
R1	1,0 mL	1,0 mL	1,0 mL	1,0 mL	1,0 mL
Homogenize and incubate the test tubes for 1 minute at 37 °C					
R2	250 µL	---	250 µL	---	250 µL

- Homogenize and incubate at 37 °C for 4 minutes.
- Measure the Blank's, Sample Blank's, Sample's and Standard's absorbance against the purified water at 550 nm (505-580 nm). The final reaction is stable for 30 minutes.

B) CALCULATIONS

Total Bilirubin (mg/dL) = $\frac{\text{Sample's } \Delta A - \text{Blank}}{\text{Calibrator's } \Delta A - \text{Blank}} \times \text{Calibrator Concentration (mg/dL)}$

In which:

ΔA = Sample's or Calibrator's Absorbance – Sample or Calibrator Blank's Absorbance

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{Calibrator Concentration (mg/dL)}}{\text{Calibrator's } \Delta A - \text{Blank}}$

Total Bilirubin (mg/dL) = (Sample's ΔA – Blank) x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

Bilirubin is conjugated in the liver and excreted in bile, resulting from the breakdown of hemoglobin in the reticuloendothelial system. Its dosage assists in the diagnosis of hepatobiliary and hemolytic diseases. There are inherited or acquired diseases that are responsible for altering the stages involved in its production, uptake, storage, metabolism and excretion, changing the balance between unconjugated (indirect), conjugated (direct) or both (total) bilirubin, causing hyperbilirubinemia. In neonates, the most common form of increased indirect bilirubin is in physiological jaundice, which may trigger encephalopathy (kernicterus). Direct bilirubin increases in the body due to inherited liver diseases such as Dubin-Johnson and Rotor diseases, damage to hepatocytes (viral, toxic or alcoholic), biliary obstruction (lithiasis or neoplasms), acute or chronic hepatitis and toxic reactions to various medications (such as chlorpromazine). Indirect bilirubin increases with hemolytic anemia, autoimmune hemolysis, blood transfusion, hematoma reabsorption, ineffective erythropoiesis, inherited diseases such as Gilbert and Crigler-Najar diseases and neonatal jaundice.

PERFORMANCE CHARACTERISTICS

Operating range
0.12 to 50.00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.07 mg/dL	0.12 mg/dL

Analytical Specificity		
Hemoglobin	Triglycerides	Ascorbic Acid
200 mg/dL	1000 mg/dL	20 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy

Number of Samples	40 in duplicate
Regression Equation	y = 0,9625x + 0,0229
Correlation Coefficient (R)	0,9983

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 1.40 mg/dL, 2.50 mg/dL and 20.00 mg/dL was -2.10%, -2.80% and -3.60%, respectively.

Precision:

Within-run precision was determined with 5 replicates for 10 days; between-run precision was determined in 3 condition, with 3 replicates for 5 days.

Samples (mg/dL)	Within-Run Precision		Samples (mg/dL)	Between-Run Precision	
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)
0.91	0.02	1.80	0.91	0.02	1.80
6.02	0.10	1.70	5.98	0.12	2.00
10.80	0.10	0.90	10.80	0.08	0.80

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Newborns: 1 day	Up to 5.1 mg/dL
Newborns: 1 to 2 days	Up to 7.2 mg/dL
Newborns: 3 to 5 days	Up to 10.3 mg/dL
Children, Adolescents and Adults	Up to 1.2 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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