CKNAC VET

CKNAC VET/CKNAC VET
Ref. 90.024.00

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ENGLISH

INTENDED USE




Kit intended to determine the activity of total Creatine Kinase in serum samples. Veterinary Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 14 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

WORKING PRINCIPLE

Method: Kinetic UV
Creatine Kinase (CK) catalyzes the phosphorylation of ADP by creatine-phosphate to form creatine and ATP. The enzymatic activity is determined by coupling the reactions of Hexokinase (HK) and Glucose-6-phosphate Dehydrogenase (G6P-DH) and is proportional to the rate of formation of NADPH, which is spectrophotometrically measured at 340 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Serum	4 to 8 °C	1 month
	20 a 25 °C	4 hours
	-20 °C	1 month

PRODUCT DESCRIPTION

R 1

Imidazole buffer ≥ 50 mmol/L, NADP ≥ 1 mmol/L, D-glicose ≥ 20 mmol/L, N-Acetyl Cysteine ≥ 20 mmol/L, hexokinase ≥ 2000 U/L, activators, preservative.

X

R 2

Creatine phosphate ≥ 100 mmol/L, AMP ≥ 5 mmol/L, Di-adenosine-5-pentafosphate ≥ 20 μmol/L, ADP ≥ 5 mmol/L, glicose-6-phosphate dehydrogenase ≥ 4000 U/L, activators, preservatives.

X

The determination of CKNAC activity is traceable to the reference method IFCC

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal VET	90.039.00
Normal Control – Quantinorm VET	90.040.00
Pathological Control – Quantialt VET	90.041.00

REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

- Heat the work reagent at 37 °C for 3 minutes.
- Zero the reading equipment with purified water filter at 340 nm.
- Pipette in the assay tube:

	Sample
Sample	40 μL
Work Reagent	1.0 mL

- Homogenize and insert it immediately in the thermostated cuvette at 37 °C.
- After 2 minutes, measure the initial absorbance A₀ and read it again every minute, for 3 minutes (A₁, A₂ and A₃, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

The CKNAC activity is calculate by multiplying the ΔA/min by the factor below:

$$\text{CKNAC (U/L)} = \Delta A/\text{min} \times 4127$$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaitda.ind.br

Note: to calculate the Factor for automatic analyzers use the absorption coefficient of NADPH = 6.3.

C) INTERPRETATION

CK is a dimeric enzyme that occurs in four different forms: a mitochondrial isoenzyme and the cytosolic isoenzymes CK-MM (skeletal muscle), CK-BB (brain) and CK-MB (myocardial). In addition to muscle tissue, CK can be located, in smaller quantities, in the kidney, brain, diaphragm, gastrointestinal tract, uterus and bladder. The determination of CK activity and its isoenzymes are used in the diagnosis and monitoring of myocardial infarction and myopathies such as Duchenne muscular dystrophy. In neurogenic muscular diseases there is no change of serum enzyme activity. It is very high in malignant hyperthermia.

PERFORMANCE CHARACTERISTICS

Operating range	
8,35 to 1000 U/L	
Sensitivity	
Detection Limit	Quantification Limit
4,43 U/L	8,35 U/L

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	200 mg/dL	10 mg/dL	1000 mg/dL
Felinos	500 mg/dL	40 mg/dL	1000 mg/dL
Equinos	100 mg/dL	40 mg/dL	1000 mg/dL
Bovinos	100 mg/dL	40 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	30 in duplicate
Regression Equation	$y = 1.0191x - 3.4607$
Correlation Coefficient (R)	0.9971
Accuracy – Felines	
Number of Samples	17 in duplicate
Regression Equation	$y = 0.9397x + 6.3471$
Correlation Coefficient (R)	0.9983
Accuracy – Equines	
Number of Samples	18 in duplicate
Regression Equation	$y = 0.981x - 0.2925$
Correlation Coefficient (R)	0.9911
Accuracy – Bovines	
Number of Samples	10 in duplicate
Regression Equation	$y = 0.9904x - 0.0654$
Correlation Coefficient (R)	0.9989

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
74.39	80	1.756	2.40	2.279	3.10
176.18	80	1.502	0.90	2.302	1.30
269.53	80	5.864	2.20	7.050	2.60

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Canines	1,15 – 28,40 U/L
Felines	7,2 – 28,2 U/L
Equines	2,4 – 23,4 U/L
Bovines	4,8 – 12,1 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

CKNAC (U/L) x 0.017 = CKNAC (µKat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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- THOMAS, L.; et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. *J Lab Med*, v. 29, p. 301-308, 2005.
- SCHUMANN, G.; et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C. Part 2. Reference Procedure for the Measurement of Catalytic Concentration of Creatine kinase. *Clin. Chem. Lab. Med.* v. 40 p. 635-642, 2002.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACCP Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. *Clinical biochemistry of domestic animals*. 5th. ed. San Diego: Academic Press, c1997 932p.

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