		<b>Responsável Técnico:</b> Dr. Gilson Serio Pizzo CRF MG – 5310 <b>Anvisa 80027310212</b>
		<b>Cálcio ASX</b> Calcium ASX / Cálcio ASX Ref. 12.002.00

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**ENGLISH**  
**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

**INTENDED USE**  
 Kit intended to determine calcium in serum, plasma and urine samples. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 15 to 30 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30°C) are followed.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**  
**Method:** Arsenazo III  
 Calcium forms a blue colored complex with Arsenazo III that can be spectrophotometrically measured at 650 - 660 nm. Color intensity is proportional to the calcium concentration in the sample.

**SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY**  
**Sample Type:** serum, heparin plasma and urine.  
**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.  
**Preservation:**

	Temperature	Sample Stability
Serum and Plasma	4 to 8 °C	14 days
	-20 °C	6 months

**Preparation:**  
**Urine:** use a sample collected within 24 hours in a bottle containing 20 to 30 mL of a 6 mol/L HCl solution to prevent the precipitation of calcium salts. Homogenize well, measure its volume and proceed with the test.

**PRODUCT DESCRIPTION**

<b>R 1</b>	MES buffer ≥ 50 mmol/L; Arsenazo III ≥ 0,1 mmol/L
	<b>STD</b>

Calcium chloride in concentration equivalent to 10 mg/dL; preservative. Traceable to the reference material NIST 3109a.

**QUALITY CONTROL**  
 The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Autocal	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 650 – 660 nm.
- Glass pipettes and/or automatic, clock or chronometer.

**TEST PROCEDURE, CALCULATION, AND INTERPRETATION**

**A) TEST PROCEDURE**  
 1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 15 to 30 °C for 2 minutes.  
 3. Measure the Sample's and Standard's absorbance against the blank at 650 – 660 nm. The final reaction is stable for 20 minutes.

**B) CALCULATIONS**  
 Cálcio (mg/dL) =  $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$   
**Calculations with the Calibration Factor:**

Calibration Factor =  $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$   
 Cálcio (mg/dL) = Sample Absorbance x Calibration Factor

**Calculations for urine samples:**  
 Cálcio urine (mg/dL) =  $\frac{\text{Sample Absorbance} \times \text{Standard Concentration (mg/dL)}}{\text{Standard Absorbance}}$   
 Cálcio (mg/24 hours) =  $\frac{\text{Cálcio (mg/dL)} \times \text{24 hours urine volume (mL)}}{100}$

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnica24.ind.br](http://www.biotechnica24.ind.br)

**C) INTERPRETATION**  
 Calcium is the most predominant cation in the organism. The skeleton contains 99% of body calcium with an unknown structure and composition similar to that of hydroxyapatite. The remainder of the calcium is distributed through the tissues and extracellular fluid where it performs specific functions in vital processes, such as: blood coagulation, neuromuscular conduction, cardiac and skeletal muscle excitability, enzymatic activation, and preservation of cell membrane permeability and integrity. Serum calcium levels are controlled by parathyroid hormone (PTH), calcitonin, and vitamin D. An imbalance in some of these modulators causes changes in serum and body calcium levels. Serum increases of PTH or vitamin D are generally associated with hypercalcemia. Increased levels of serum calcium can be observed in multiple myeloma and other neoplastic diseases. Hypocalcemia can be observed for example in hypoparathyroidism, nephrosis and pancreatitis.

**PERFORMANCE CHARACTERISTICS**

Operating range
0.94 a 20.00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.45 mg/dL	0.94 mg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	15 mg/dL	900 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 0.988x + 0.081
Correlation Coefficient (R)	0.9979

By applying the regression equation, the total systematic error estimated for decision levels of 5.0 mg/dL and 16.0 mg/dL was 0,42% and 0,69%, respectively.  
**Precision:**

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
8.21	80	0.036	0.40	0.037	0.50
12.22	80	0.054	0.40	0.054	0.40
17.44	80	0.068	0.40	0.074	0.40

CV: Coefficient of variation; SD: Standard deviation

## RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

## REFERENCE RANGES

	mg/dL	mmol/L
Serum and Plasma	8,6 - 10,3	2,15 - 2,57
	mg/24 hours	mmol/24 hours
Urine	Up to 300	Up to 7,49

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Calcium (mg/dL) x 0,25 = Calcium (mmol/L)

## WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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