Cistatina C

Cystatin C/ Cistatina C
Ref. 20.026.00

Responsável Técnico:
Dr. Gilson Serio Pizzo
CRF MG – 5310
Anvisa 80027310305

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ENGLISH

 **BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

INTENDED USE

Kit intended to determine Cystatin C in serum and plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Immunoturbidimetric

Cystatin C from sample binds to the reagent's anti-Cystatin C antibodies in an agglutination reaction. The formed antigen-antibody complex increases the turbidity of the solution, which is spectrophotometrically determined at 546 nm. The intensity of the turbidity is proportional to the Cystatin C concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma (EDTA or heparin).

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	1 month
	-20 °C	6 months

PRODUCT DESCRIPTION

R 1	Tham 0,01 mmol/L, PEG 6000 25,0 g/L, sodium chloride 16,0 g/L and preservative 1,0 g/L.	
R 2	Tham 0,01 mmol/L and latex particles coated with anti-Cystatin C antibodies 310 mL/L.	
CAL	Sodium chloride 100 mmol/L, BSA 20 g/L, sucrose 10%, preservative 0,03% and cystatin C. The calibration set consists of 5 calibration points. Consult the concentration in the product's label.	
CONTROL 1	Sodium chloride 100 mmol/L, BSA 20 g/L, sucrose 10%, preservative 0,03% and cystatin C $\approx$ 1,77 mg/L.	
CONTROL 2	Sodium chloride 100 mmol/L, BSA 20 g/L, sucrose 10%, preservative 0,03% and cystatin C $\approx$ 4,80 mg/L.	

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrator – Cystatin C CAL	REF	20.026.00
Level 1 Control – Cystatin C CONTROL 1		20.026.00
Level 2 Control – Cystatin C CONTROL 2		20.026.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 546 nm.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.ind.br

Reaction Type	Two points / Positive
Wavelength	546 nm
Temperature	37 °C
Sample volume*	3 µL
R1 Volume*	240 µL
Sample + R1 Incubation Time	300 seconds
R2 Volume*	60 µL
R1 + Sample + R2 Incubation Time	300 seconds
First Reading (A1)	Immediately
Second Reading (A2)	300 seconds
Calibration**	Calibration curve

* Volumes can be modified as long as the stablished proportion is maintained.

** The calibration set consists of 5 calibration points. Consult the concentration in the product's label.

B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

C) INTERPRETATION

Cystatin C is a cysteine-proteinase inhibitor produced by all nucleated cells. Due to its size, cystatin C is freely filtered by the renal glomerulus and is catabolized by the renal tubular epithelial cells. Therefore, it does not return to the bloodstream. Besides, cystatin C is produced at a constant rate, so it can be used as an indicator of glomerular filtration rate (GFR). As the cystatin levels are less dependent on age, sex, race and muscle mass than creatinine, this analyte is a better option to evaluate renal function.

PERFORMANCE CHARACTERISTICS

Operating range
0.20 to 8.00 mg/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Detection Limit
0.15 mg/L

Analytical Specificity			
Billirubin	Hemoglobin	Triglycerides	Rheumatoid Factor
50 mg/dL	500 mg/dL	500 mg/dL	120 IU/L

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y = 0.999x - 0.0990
Correlation Coefficient (R)	0.9980

The product's accuracy was evaluated with a cystatin C sample, at the decision level of 1,15 mg/L, with a total systematic error of 8,71%.

Precision:

Within-run precision was determined with 3 lots and 10 replicates in an analytical run; between-run precision was determined with 3 lots and 3 replicates in an analytical run.

Samples (mg/L)	Within-Run Precision		Between-Run Precision	
	SD (mg/L)	CV (%)	SD (mg/L)	CV (%)
1.57	0.055	3.50	0.052	5.58

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Serum and Plasma	0.55 - 1.15 mg/L
------------------	------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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