## Cloretos

Clorídeo / Cloruros  
Ref. 12.003.00

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#### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE ENGLISH VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

#### INTENDED USE

Kit intended to determine chloride in serum, plasma, urine and cerebrospinal fluid (CSF) samples. Diagnostic use only.


#### STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30°C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Mercury Thiocyanate

Chloride ions from sample react with mercury thiocyanate to form mercuric chloride and free thiocyanate ions. The latter react with ferric ions to form ferric thiocyanate, whose red color can be spectrophotometrically measured at 500 nm. Color intensity is proportional to the chloride concentration in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

**Sample Type:** serum, heparin plasma, urine and cerebrospinal fluid (CSF).

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Sample Stability
Serum*, Plasma* and Cerebrospinal Fluid	4 to 8 °C	7 days
	20 a 25 °C	7 days
	-20 °C	1 year
24 Hours Urine	2 to 8 °C	14 days
	15 a 30 °C	14 days
	-20 °C	2 months

\* Serum and plasma samples must be separated within 1 hour after collection.

#### Preparation urine 24 hours:

**Urine:** use a sample collected within 24 hours, with no need to add preservatives. Homogenize and measure its volume. Centrifuge an aliquot of the urine and collect the supernatant. Dilute an aliquot of the supernatant in a 1:2 ratio with purified water and proceed with the assay. Multiply the obtained result by 2. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

**Cerebrospinal fluid:** centrifuge for 10 minutes at 3000 rpm and use the supernatant in the assay.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Mercury thiocyanate $\geq 1.0$ mmol/L, ferric nitrate $\geq 20$ mmol/L, nitric acid $\geq 0.1\%$ v/v, sulfuric acid $\geq 0.1\%$ v/v, mercury nitrate $\geq 0.001\%$ v/v.	
<b>STD</b>	Chloride ions at a concentration equivalent to 100 mEq/L; preservative. Traceable to the certified reference materials NIST SRM 3182 and NIST SRM 999b.	

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Autocal H	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

**REF**

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 500 nm (420 – 520 nm).
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	5 µL	-
Sample	-	-	5 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 15 to 30 °C for 5 minutes.

3. Measure the Sample's and Standard's absorbance against the blank at 500 nm (420 – 520 nm). The final reaction is stable for 30 minutes.

##### B) CALCULATIONS

Chloride (mEq/L) =  $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mEq/L)}$

**Calculations with the Calibration Factor:**

Calibration Factor =  $\frac{\text{STD Concentration (mEq/L)}}{\text{STD Absorbance}}$

Chloride (mEq/L) = Sample Absorbance x Calibration Factor

**Calculations for urine samples:**

Chloride (mEq/24 hours) = Chloride (mEq/L) \* x Urinary Volume 24 hours (L)

\*The supernatant of the urine sample must be diluted in a 1:2 ratio. To determine the chloride concentration in the urine sample (mEq/L), multiply the obtained result by 2.

**Automação:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br)

##### C) INTERPRETATION

Chloride is the major extracellular anion. Together, sodium and chloride represent the bulk of the osmotically active constituents of plasma. The chloride ions from the feed are almost completely absorbed in the intestinal tract. It is significantly involved in: maintenance of water distribution, osmotic pressure and anion-cation balance in the extracellular fluid (ECF). Determination of plasma chloride concentration is useful in the differential diagnosis of acid-base disorders and is essential for calculating the anionic range. Chloride fluctuations in serum or plasma have little clinical consequence, but serve as signs of a basic disorder in fluid homeostasis and acid-base balance, and may assist in differentiating the cause of these disorders. Hyperchloremia occurs in dehydration, renal tubular acidosis (RTA), acute renal failure, metabolic acidosis associated with prolonged diarrhea and loss of sodium bicarbonate, diabetes insipidus, conditions of adrenocortical hyperfunction, overload of saline solutions and respiratory alkalosis as a result of renal excretion of HCO<sub>3</sub><sup>-</sup>. Hypochloremia is observed in metabolic acidosis caused by increased production or decreased excretion of organic acids (eg, diabetic ketoacidosis and renal insufficiency), persistent gastric secretion, prolonged vomiting, hypochloremic alkalosis, and total body chloride depletion with corresponding retention of HCO<sub>3</sub><sup>-</sup>. Other associated conditions include salt loss nephritis,

Addison's disease, hyperaldosteronism, and bromide poisoning. The excess chloride is excreted in the urine and is also eliminated in sweat. Usually urinary excretion of chlorine approaches dietary intake. Physiologically increased levels of post-menstrual diuresis and decreased pre-menstrual water retention of salt and water are observed.

#### PERFORMANCE CHARACTERISTICS

Operating range	
Serum / Plasma / Urine	66.87 to 150.00 mEq/L
Cerebrospinal Fluid	79.12 a 132.90 mEq/L

For concentrations above the operating range, dilute the sample with purified water, proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity		
Sample Type	Detection Limit	Quantification Limit
Serum and Plasma	3.05 mEq/L	66.87 mEq/L
Cerebrospinal Fluid	-	79.12 mEq/L

Analytical Specificity			
Sample Type	Hemoglobin	Bilirubin	Triglycerides
Serum and Plasma	500 mg/dL	40 mg/dL	1250 mg/dL
Cerebrospinal Fluid	500 mg/dL	-	-

Concentrations of interfering substances up to the values presented above do not cause significant changes in the results. Any concentration of membrane stabilizer (Transfix\*) interferes with the determination. For drugs, consult the recommended bibliographic reference (Young, 2000).

Sample Type	Accuracy	
	Number of Samples	40 in duplicate
Serum and Plasma	Regression Equation	$y = 0.987x + 1.268$
	Correlation Coefficient (R)	0.9972
	Number of Samples	20 in duplicate
Cerebrospinal Fluid	Regression Equation	$Y = 1.0103x + 5.8554$
	Correlation Coefficient (R)	0.9952
	Number of Samples	18 in duplicate
Urine	Regression Equation	$Y = 0.997x + 0.477$
	Correlation Coefficient (R)	0.9995

**Serum and Plasma:** by applying the regression equation, the total systematic error estimated for decision levels of 70.0 mEq/L and 120.0 mEq/L, was 0.51% and -0.24%, respectively. **Cerebrospinal Fluid:** by applying the regression equation, the total systematic error estimated for decision levels of 80.0 mEq/L, 107.7 mEq/L and 130.0 mEq/L, was 8.35%, 6.46% and 5.53%, respectively.

**Urine:** by applying the regression equation, the total systematic error estimated for decision levels of 120 mEq/24h e 220,0 mEq/24h was 0,0,9% and 0,08%, respectively.

**Precision:**

Sample Type	Samples (mEq/L)	Replicates	Within-Run Precision		Total Precision	
			SD (mEq/L)	CV (%)	SD (mEq/L)	CV (%)
Serum and Plasma	86.14	80	0.769	0.90	0.836	1.00
	95.21	80	0.804	0.80	0.852	0.90
	135.92	80	0.815	0.60	0.893	0.70
Cerebrospinal Fluid	88.51	20/40	1.08	1.22	1.21	1.36
	99.32	20/40	3.26	2.10	2.51	2.51
	125.83	20/40	2.65	2.10	2.19	1.74

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- R1 is toxic because it contains Mercury thiocyanate.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- R1 is toxic because it contains Mercury thiocyanate.

#### REFERENCE RANGES

Serum and Plasma	98 - 107 mEq/L
24 Hours Urine	110 - 250 mEq/24 hours
Cerebrospinal Fluid	112 - 123 mEq/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Chloride (mEq/L) x 1.0 = Chloride (mmol/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

• All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

• To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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