## Cobre Plus

Copper Plus / Cobre Plus  
Ref. 12.013.00

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### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

### INTENDED USE

Kit intended to determine copper ions in serum samples. Diagnostic use only.

### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing. Do not freeze.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed. Store protected from light.
- Do not use reagents whose shelf life has expired.

### WORKING PRINCIPLE

Around 95% of body's copper ions is found bounded to ceruloplasmin. In acidic pH, ions copper dissociate from ceruloplasmin and react with the chromogen di-Br-PAESA, forming a colored complex that can be measured at 580 nm. Color intensity is proportional to the concentration of ions copper in sample.

### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum.

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

### Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	2 weeks
	-20 °C	1 year

### PRODUCT DESCRIPTION

<b>R 1</b>	Acetate buffer < 500 mM; 3,5 Di-Br-PAESA pH > 3,0; surfactants and stabilizers.	<b>X</b>
<b>STD</b>	Copper solution and stabilizers. Concentration on the bottle label. Value attributed according to the Atomic Absorption methodology.	<b>X</b>
<b>CONTROL 1</b>	Copper solution and stabilizers. Value range on bottle label.	<b>X</b>
<b>CONTROL 2</b>	Copper solution and stabilizers. Value range on bottle label.	<b>X</b>

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the control below:

Quantialt Control	<b>REF</b>	13.004.00
Cobre Plus CONTROL 1	<b>REF</b>	12.013.00
Cobre Plus CONTROL 2	<b>REF</b>	12.013.00

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 580 nm (560 - 630 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Before starting, let the reagents on the bench for at least 30 minutes until they reach room temperature.
2. Pipette in the test tubes:

	Blank	STD	Sample
R1	1000 µL	1000 µL	1000 µL
Incubate for 3 minutes at 37°C and then add:			
STD	---	50 µL	---
Sample	---	---	50 µL

3. Homogenize and incubate the test tubes at 37 °C for 6 minutes.
4. Measure the absorbance of the STD and sample against the reagent blank.

### B) CALCULATIONS

Copper (µg/dL) =  $\frac{\text{Sample absorbance} \times \text{Standard Concentration}}{\text{Standard absorbance}}$  (µg/dL)

#### Example:

Standard Concentration = 200 µg/dL  
Absorbance of the sample = 0,161  
Absorbance of standard = 0,230

Copper (µg/dL) =  $\frac{0,161 \times 200}{0,230}$  = 140,0 µg/dL

#### Calculations with the Calibration Factor:

Calibration factor =  $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$

Copper (µg/dL) = Sample Absorbance x Calibration factor

#### Example:

Calibration factor =  $\frac{200}{0,230}$  = 869,56

Copper (µg/dL) = 0,161 x 869,56 = 140,0 µg/dL

**Note:** Standard concentration on the bottle label.

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

### C) INTERPRETATION

Ion's copper are cofactor to various enzymatic reactions, besides being part of iron transport and absorption metabolism. Reduced levels of copper are correlated with diminished serum proteins, which is observed in the following cases: malnutrition, incorrect absorption (celiac disease), loss of protein through feces or urine (nephrotic syndrome) or Wilson's disease. Increased levels of copper can be observed during pregnancy, chronic or acute infections, myocardial infarction, hyperthyroidism and hematological diseases.

### PERFORMANCE CHARACTERISTICS

Operating range
6,0 a 500,0 µg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0,65 µg/dL	5,50 µg/dL

Analytical Specificity		
Bilirubin	Hemolysis	Lipemia
40 mg/dL	Interferes with the dosage	

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	50 in duplicate
Regression Equation	$y = 0,9627x + 5,1975$
Correlation Coefficient (R)	0,9978

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 80,0 µg/dL, 180,0 µg/dL and 270,0 µg/dL was 2,8%, -0,8% and -1,8%, respectively.

**Precisión:**  
Within-run precision was determined with 5 replicates for 10 days; between-run precision was determined in 3 conditions, with 3 replicates for 5 days.

Samples (µg/dL)	Within-Run Precision		Samples (µg/dL)	Between-Run Precision	
	SD (µg/dL)	CV (%)		SD (µg/dL)	CV (%)
16.24	0.39	2.4%	16.33	0.31	1.9%
126.72	1.04	0.8%	126.89	1.14	0.9%
245.99	2.59	1.1%	246.73	1.73	0.7%

CV: Coefficient of variation; SD: Standard deviation.

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination.
- O R1 tends to precipitate at low temperatures, therefore wait for the reagent to reach room temperature before use, as specified in the test procedure.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Male	70 – 140 µg/dL
Female	80 – 155 µg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646.

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br)
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br)

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