Colesterol

Cholesterol / Colesterol
Ref. 10.004.00 / FD38.048.00 / SD39.057.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine cholesterol in serum samples. Diagnostic use only.



STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light.
- Reagent ready for use. Due to STD's volatility, the flask must be immediately closed after its use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- The Dedicated and Semi Dedicated Reagents presentations have onboard stability of 28 days.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic - Trinder

Cholesterol esters from sample are hydrolyzed by the enzyme cholesterol esterase to form free cholesterol. The enzyme cholesterol oxidase, in the presence of oxygen, catalyzes the oxidation of the free cholesterol to form hydrogen peroxide. In the presence of the latter and of 4-aminoantipyrine, peroxidase (POD) catalyzes the oxidation of phenol to form a pink compound (4-(p-benzoquinon-monoimino)-fenazon) that can be spectrophotometrically measured at 505 nm. Color intensity is proportional to the cholesterol concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	7 days
	-20 °C	3 months

PRODUCT DESCRIPTION

Pipes buffer ≥ 20 mmol/L; 4-aminoantipyrine ≥ 0,1 mmol/L; phenol ≥ 1,0 mmol/L; peroxidase ≥ 1000 U/L; cholesterol oxidase ≥ 200U/L; cholesterol esterase ≥ 200 U/L; lipoprotein lipase ≥ 200U/L; activators; detergents, stabilizers and preservative.

R 1

STD

Butyl alcohol ≥ 5% v/v; detergent, stabilizers and Cholesterol in concentration equivalent to 200 mg/dL. Traceable to reference material NIST 1951b.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator – Autocal	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 505 nm (490 – 510 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 10 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 505 nm (490 – 510 nm). The final reaction is stable for 25 minutes.

B) CALCULATIONS

Colesterol (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Colesterol (mg/dL) = Sample Absorbance x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicalda.ind.br

C) INTERPRETATION

Cholesterol is a steroid with a secondary hydroxyl group on the C₃ position and it is an essential component of the membrane of all cells. Approximately 75% of it is synthesized endogenously, mainly in the liver and in the intestinal wall, and the remaining is obtained through the diet. Most peripheral cells depend on the exogenous distribution of cholesterol by lipoproteins. Its elevation can be associated with atherosclerotic risk or other metabolic disorders of lipids and lipoproteins.

PERFORMANCE CHARACTERISTICS

Operating range
4,60 to 800,00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1,74 mg/dL	4,60 mg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	15 mg/dL	1250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 0.992x - 4.48
Correlation Coefficient (R)	0.9908

By applying the regression equation, the total systematic error estimated for decision levels of 200 mg/dL and 500 mg/dL was -3.04% and -1.69%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
135,77	80	1,20	0,90	1,36	1,00
295,52	80	2,45	0,80	2,98	1,00
595,42	80	2,62	0,40	2,86	0,50

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

According to the V Brazilian directive Cholesterol levels can be evaluated based on age as follows.

Age	Total Cholesterol Values (mg/dL)	Interpretation
2 to 19 years	< 150	Desirable
	150 – 169	Boarder Range
	≥ 170	High
20 years or more	< 200	Desirable
	200 – 239	Boarder Range
	> 240	High

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Cholesterol (mg/dL) x 0,026 = Cholesterol (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

• To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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