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Dr. Gilson Serio Pizzo

CRF MG – 5310

Ref. 20.006.00

Anvisa 80027310173

Complemento C4

C4 Complement / Complemento C4

Ref. 20.006.00

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**ENGLISH**

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

**INTENDED USE**

Kit intended to determine C4 complement (C4) in serum and plasma samples. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**

**Method:** Immunoturbidimetric.

Anti-complement C4 antibodies form insoluble complexes when mixed with samples containing C4. The turbidity intensity is proportional to the C4 concentration of the sample.

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**

**Sample Type:** serum and plasma (heparin or EDTA).

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

|                | Temperature         | Stability Period   |
|----------------|---------------------|--------------------|
| Serum / Plasma | 4 to 8 °C<br>-20 °C | 8 days<br>3 months |

**PRODUCT DESCRIPTION**

R1

Goat anti-complement C4 antibodies; phosphate buffer; preservative.

**QUALITY CONTROL**

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and control below:

|                             |           |
|-----------------------------|-----------|
| Multicalibrator Plus        | 21.011.00 |
| Medium Multicalibrator Plus | 21.003.00 |

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 340 – 405 nm.
- Test Tubes, glass pipettes and/or automatic.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

**Automation:** This product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnica Ltda.ind.br](http://www.biotecnica Ltda.ind.br)

| Reaction Type             | End Point         |
|---------------------------|-------------------|
| Wave length               | 340 nm            |
| Temperature               | 37 °C             |
| Sample Volume *           | 5,5 µL            |
| R1 Volume *               | 220 µL            |
| R1 Incubation Time        | 180 seconds       |
| R1 + Sample Reaction Time | 120 seconds       |
| Calibration **            | Calibration curve |

\* Sample and R1 volumes can be modified as long as the established ratio is maintained.

**\*\* CALIBRATION CURVE:**

Dilute the Multicalibrator Plus in a 9,0 g/L NaCl solution as follows:

| Dilution          | 1   | 2   | 3    | 4   | 5    | 6   |
|-------------------|-----|-----|------|-----|------|-----|
| Calibrator (µL)   | -   | 10  | 25   | 50  | 75   | 100 |
| NaCl 9,0 g/L (µL) | 100 | 90  | 75   | 50  | 25   | -   |
| Dilution Factor   | 0   | 0.1 | 0.25 | 0.5 | 0.75 | 1   |

Multiply the C4 complement concentration indicated in the Multicalibrator Plus instruction of use by the corresponding dilution factor to obtain the C4 concentration of each dilution. The C4 concentration of the Multicalibrator Plus is traceable to the Standard Reference Material ERM DA470k / IFCC.

**B) CALCULATIONS**

The equipment automatically calculates the analyte concentration in each sample.

**C) INTERPRETATION**

C4 is a component of the essential complement for pathway activation of this system. It is a protein synthesized by the liver and to a lesser extent by monocytes and other tissues. Increases and decreases in concentration are of great clinical significance. Increased levels of C4 concentration are related to acute phase responses (trauma, inflammation). Decreases in concentration are related to deficiencies of genetic origin (autoimmune or collagen vascular diseases, particularly Systemic Lupus Erythematosus), or acquired deficiency as a consequence of the formation of immune complexes, autoimmune hemolytic anemia, and autoimmune.

**PERFORMANCE CHARACTERISTICS**

| Operating range   |
|-------------------|
| 3.05 to 107 mg/dL |

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

| Sensitivity          |
|----------------------|
| Quantification Limit |
| 3.05 mg/dL           |

| Analytical Specificity |           |               |
|------------------------|-----------|---------------|
| Hemoglobin             | Bilirubin | Triglycerides |
| 75 mg/dL               | 7.5 mg/dL | 750 mg/dL     |

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy                    |                      |
|-----------------------------|----------------------|
| Number of Samples           | 40 in duplicate      |
| Regression Equation         | y = 0.9952x + 0.1418 |
| Correlation Coefficient (R) | 0.9864               |

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 10 mg/dL, 25 mg/dL and 65 mg/dL was 0.9%, 0.1% and -0.3%, respectively.

**Prozone High Dose Effect:**

Prozone effect does not occur up to 239.14 mg/dL of C4 complement.

**Precision:**

The intra-assay precision studies were performed in 5 replicates for 10 days; the inter-assay studies were performed in 3 conditions, with 3 replicates for 5 days.

| Samples (mg/dL) | Intra-assay Precision |        | Samples (mg/dL) | Inter-assay Precision |        |
|-----------------|-----------------------|--------|-----------------|-----------------------|--------|
|                 | SD (mg/dL)            | CV (%) |                 | SD (mg/dL)            | CV (%) |
| 19.59           | 0.74                  | 3.8    | 18.90           | 0.92                  | 4.8    |
| 34.09           | 0.83                  | 2.4    | 33.56           | 0.91                  | 2.7    |
| 51.30           | 1.06                  | 2.1    | 52.15           | 1.27                  | 2.4    |

CV: Coefficient of variation; SD: Standard deviation

**RESIDUAL RISKS, WARNINGS AND PRECAUTIONS**

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

**REFERENCE RANGES**

|                      |               |
|----------------------|---------------|
| Adults (20-60 years) | 10 a 40 mg/dL |
|----------------------|---------------|

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

BIOTÉCNICA IND.COM. LTDA - CNPJ: 025340690001/20, Avenida Washington Ribeiro, 200 - Industrial Miguel de Luca  
CEP 37072-030 Varginha-MG BRASIL Tel/fax: +55 35 3214 4646 [www.biotecnica.ind.br](http://www.biotecnica.ind.br)

Revisão/Revision/Revisión: 02 – 03/04/2024

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) Or calling for +55 35 3214 4646.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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- PRICE, C. P.; et al. Ann Clin Biochem. v. 20, p. 1-14, 1983.
- DATI, F. et al; Eur J Clin Chem Clin Biochem. v. 34, p. 517-520, 1996.
- TIEZ, Textbook of Clinical Chemistry. 3 ed. Burtis CA, Ashwood ER. WB Saunders Co, 1999.
- FRIEDMAN AND YOUNG, Effects of the disease on clinical laboratory tests, 3 ed., AACC Press, 1997.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.BURTIS, C.A.; ASHWOOD, E.R.; BRUNS, D.E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.

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