## Creatinine VET

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#### ENGLISH

#### INTENDED USE

Kit intended to determine creatinine in serum and plasma samples. Veterinary diagnostic use only.

#### STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Work Reagent:** mix in the proportion of 1 part of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 21 days, as long as the preparation instructions and the recommended storage conditions are followed and stored in a closed amber plastic bottle, protected from light (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

#### WORKING PRINCIPLE

**Method:** Picrate

The creatinine in the sample reacts with the Picrate in alkaline medium to form an orange-red complex that can be spectrophotometrically measured at 500 nm. Color intensity is proportional to the creatinine concentration in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

**Sample Type:** serum, plasma (EDTA and heparin).

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum	4 to 8 °C	7 days
	-20 °C	3 months

#### PRODUCT DESCRIPTION

<b>R 1</b>	Picric acid ≥ 10 mmol/L.	
<b>R 2</b>	Sodium Hydroxide ≥ 200 mmol/L, amino acetic acid ≥ 10 mmol/L, surfactant.	
<b>STD</b>	Stabilizer; preservative; Creatinine in concentration equivalent to 2 mg/dL. Traceable to the reference material NIST 914a.	

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal Vet	90.039.00
Normal Control – Quantinorm Vet	90.040.00
Pathological Control – Quantialt Vet	90.041.00

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 500 nm (490 – 510 nm).

- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

- Incubate the work reagent at 37 °C for 3 minutes.
- Zero the equipment at 500 nm with purified water.
- Pipette in the test tubes:

	STD	Sample
STD	100 µL	-
Sample	-	100 µL
Work Reagent	1,0 mL	1,0 mL

- Homogenize and insert into the thermostated cuvette at 37 °C. Start the chronometer.
- Measure the Sample's and Standard's absorbance at 30 seconds (A1) and 90 seconds (A2).

##### B) CALCULATIONS

Creatinine (mg/dL) =  $\frac{\text{Sample's (A2 - A1)} \times \text{STD Concentration (mg/dL)}}{\text{STD's (A2 - A1)}}$

##### Calculations with the Calibration Factor:

Calibration Factor =  $\frac{\text{STD Concentration (mg / dL)}}{\text{STD's (A2 - A1)}}$

Creatinine (mg / dL) = Sample's (A2 - A1) x Calibration Factor

**Automation:** this product is compatible with most types of automatic analyzers. Instrument settings are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

#### C) INTERPRETATION

Creatinine is the end product of the breakdown of phosphocreatine. It is produced endogenously and is released into body fluids at a constant rate. Its plasma concentration is maintained within narrow limits, predominantly by glomerular filtration. Consequently, both plasma creatinine concentration and renal clearance (creatinine clearance) have been used as markers of glomerular filtration rate.

#### PERFORMANCE CHARACTERISTICS

Operating range
0.21 to 10.83 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
0.20 mg/dL

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglicérides
Canines	250 mg/dL	10 mg/dL	500 mg/dL
Felines	200 mg/dL	18 mg/dL	500 mg/dL
Equines	300 mg/dL	10 mg/dL	250 mg/dL
Bovines	100 mg/dL	10 mg/dL	250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	40 in duplicate
Regression Equation	y = 0.9858x + 0.0215
Correlation Coefficient (R)	0.9987
Accuracy – Felines	
Number of Samples	20 in duplicate
Regression Equation	y = 0.998x - 0.0028
Correlation Coefficient (R)	0.9980
Accuracy – Equines	
Number of Samples	11 in duplicate
Regression Equation	y = 1.0371x - 0.0435
Correlation Coefficient (R)	0.9991
Accuracy – Bovines	
Number of Samples	10 in duplicate
Regression Equation	y = 0.9546x + 0.0565
Correlation Coefficient (R)	0.9933

#### Precision:

Determined with two runs and 20 replicates.

Samples (U/L)	Number of samples	Within-Run Precision		Samples (U/L)	Number of samples	Total Precision	
		SD (U/L)	CV (%)			SD (U/L)	CV (%)
1.08	20	0.03	3.00	1.06	40	0.04	3.80
6.99	20	0.06	0.80	6.98	40	0.07	0.90

CV: Coefficient of variation; SD: Standard deviation

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

### REFERENCE RANGES

Canine	0.5 – 1.5 mg/dL
Feline	0.8 – 1.8 mg/dL
Equine	1.0 – 2.0 mg/dL
Bovine	1.2 – 1.9 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µmol/L):

Creatinine (mg/dL) x 88,4 = Creatinine (µmol/L)

### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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