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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine creatinine in serum, plasma (EDTA and heparin) and urine samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 ° C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use. It can be observed a color alteration in R2 to orange, however it does not affect the reagent's capacity as long as it is monthly calibrated.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- The Dedicated and Semi Dedicated Reagents presentations have *onboard* stability of 28 days.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic-PAP
Creatinine is hydrolyzed by creatinine amidohydrolase, forming creatine. This suffers the action of creatine amidohydrolase, generating sarcosine and urea. Through oxidative demethylation caused by sarcosine oxidase, sarcosine is transformed into glycine, formaldehyde, and hydrogen peroxide. DHBS and 4-amino antipyrine react with hydrogen peroxide, in the presence of peroxidase, producing quinoneimine, which is spectrophotometrically determined at 510 nm. The color intensity is proportional to the concentration of creatinine in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum, plasma (EDTA and heparine) and urine.
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preparation:

	Temperature	Stability Period
Serum/Plasma	4 to 8 ° C	7 days
	-20 ° C	3 months
24 Hours Urine	4 to 8 ° C	6 days
	-20 ° C	6 months

Preparation:

Urine: work with samples collected in the period of 24 hours. Homogenize the urine, measure its volume and separate 20 mL. Adjust the pH between 7 and 9 with a 5% NaOH solution. In case of turbidity, incubate the sample for 10 minutes at 56 °C. Perform a 1:20 dilution of the aliquot with purified water and perform the assay. Multiply the obtained result by 20. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

PRODUCT DESCRIPTION

R 1	Buffer pH 7,5, preservative, stabilizers, detergent, DHBS <10 mmol/L, ascorbate oxidase <10 KU/L, catalase <500 KU/L, sarcosine oxidase <50 KU/L, creatine amidinohydrolase <50 KU/L.	R 2	Buffer pH 8,0, preservative, stabilizer, 4-aminoantipyrine <10 mmol/L, peroxidase <80 KU/L, creatinine amidohydrolase <50 KU/L.	 
------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------	---------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls:
Calibrator - Autocal H 13.002.00
Normal Control - Quantinorm 13.003.00
Pathological Controle - Quantialt 13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer for reading at 510 nm (480 - 546 nm) and 670 nm (630 - 690 nm).
- Test tubes, glass pipettes and/or automatic.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

This procedure is compatible to most types of automatic analyzers. The instrument settings are available at www.biotechnicaltda.ind.br

Reaction type	End point / Crescent
Wavelength	510 nm / 670 nm
Temperature	37 °C
Volume sample*	5 µL
Volume R1*	200 µL
Read 1 (A ₁)	306 seconds
Volume R2*	50 µL
Read 2 (A ₂)	270 seconds
Calibration	2 points

* The volumes can be modified as long as the proportion established is maintained.

B) CALCULATIONS

Calculation for Serum/Plasma:
The equipment automatically calculates the concentration of the analyte in each sample.
Calculation for Urine:
Creatinine (mg/24 horas) = $\frac{\text{Creatinine (mg/dL)} \times \text{Urinary volume (mL)}}{100}$

*The urine sample must be diluted in a 1:20 proportion. To determine the creatinine concentration in the urine sample (mg/dL), multiply the obtained result by 20.

Creatinine in mg/kg/24 hours:
Creatinine (mg/kg/24 horas) = $\frac{\text{Creatinine (mg / 24 horas)}}{\text{Weight (kg)}}$

Endogenous Creatinine Clearance:
Debug (mL/min) = $\frac{\text{Urine Creatinine (mg/dL)} \times \text{VM}^{**}}{\text{Serum creatinine (mg/dL)}}$

*The urine sample must be diluted in a 1:20 proportion. To determine the creatinine concentration in the urine sample (mg/dL), multiply the obtained result by 20.
**VM = minute volume (24 hours urine volume in mL divided by 1440 minutes).
Attention: The clearance should be corrected to the patient's body surface.

Calculation of the Body Surface (SC) without using the Nomogram:
SC (m²) = Weight (kg)^{0,425} x Heighy (cm)^{0,725} x 0,007184
Clearance:
Clearance (mL/min) = $\frac{\text{Debug (mL/min)} \times 1,73}{\text{Body Surface}}$

C) INTERPRETATION

Creatinine is the product of phosphocreatine decomposition. It is produced endogenously and released into body fluids at a constant rate and its plasma concentration is maintained within narrow limits, predominantly by glomerular filtration. Consequently, both the

plasma creatinine concentration and its renal clearance (creatinine clearance) have been used as markers of glomerular filtration rate.

PERFORMANCE CHARACTERISTICS

Operating range	
0.20 to 80.00 mg/dL	

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

	Sensitivity	
	Detection Limit	Quantification Limit
Serum	0.021 mg/dL	0.174 mg/dL
Urine	0.021 mg/dL	0.030 mg/dL

	Analytical Specificity				
	Hemoglobin	Bilirubin	Triglycerides	Creatine	Ascorbic Acid
Serum	500 mg/dL	10 mg/dL	1000 mg/dL	40 mg/dL	100 mg/dL
Urine	500 mg/dL	5 mg/dL	-	40 mg/dL	100 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine
Number of samples	45 in duplicate	19 in duplicate
Regression equation	$y = 0.9434x + 0.0259$	$y = 0.9445x + 3.8567$
Correlation coefficient (R)	0.9999	0.9999

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 0.6 mg/dL, 1.6 mg/dL and 6.0 mg/dL was -1.3%, -4.0% and 5.2%, respectively.

Urine: by applying the regression equation, the total systematic error estimated for decision levels 750 mg/24h, 1500 mg/24h and 3000 mg/24h was -5.0%, -5.3% and -5.4%, respectively.

Precision:

The intra-assay precision studies were performed in 5 replicates for 10 days, while the inter-assay studies were performed in 3 conditions, with 3 replicates for 5 days.

Serum Samples (mg/dL)	Intra-assay Precision		Serum Samples (mg/dL)	Inter-assay Precision	
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)
0.68	0.03	4.1	0.72	0.02	2.3
3.79	0.12	3.1	3.92	0.07	1.9
8.61	0.23	2.7	8.86	0.23	2.6

Urine Samples (mg/dL)	Intra-assay Precision		Urine Samples (mg/dL)	Inter-assay Precision	
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)
76.80	2.39	3.1	79.88	1.80	2.2
188.13	5.63	3.0	194.12	4.75	2.4
358.65	14.76	4.1	367.09	16.33	4.4

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Serum/Plasma (mg/dL)			
Newborn	0.31 - 0.92	9 - 11 years	0.32 - 0.64
2 week - 3 years	0.16 - 0.39	11 - 13 years	0.42 - 0.71
3 - 5 years	0.26 - 0.42	13 - 15 years	0.46 - 0.81
5 - 7 years	0.29 - 0.48	Adult (man)	0.73 - 1.18
7 - 9 years	0.34 - 0.55	Adult (woman)	0.55 - 1.02

Urine (mg/kg/24h)	Creatinine clearance (mL/min/1,73m ²)
2 - 3 years	6 - 22
>3 years	12 - 30
Adult (man)	21 - 26
Adult (woman)	16 - 22

These values are intended for orientation only. It's recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotechnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotechnica

instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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