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### ENGLISH

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

### INTENDED USE

Kit intended to determine creatinine in serum, plasma and urine samples. Diagnostic use only.

### STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- **Work Reagent:** mix in the proportion of 1 part of R1 + 1 part of R2. Homogenize gently. The working reagent is stable for 21 days, as long as the preparation conditions are followed and stored in a closed amber plastic bottle, protected from light (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

### WORKING PRINCIPLE

**Method:** Picrate  
The creatinine in the sample reacts with the Picrate in alkaline medium to form an orange-red complex that can be spectrophotometrically measured at 500 nm. Color intensity is proportional to the creatinine concentration in the sample.

### SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

**Sample Type:** serum, plasma (EDTA and heparin) and urine.  
**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

|                | Temperature | Stability Period |
|----------------|-------------|------------------|
| Serum          | 4 to 8 °C   | 7 days           |
|                | -20 °C      | 3 months         |
|                | 20 to 25 °C | 7 days           |
| 24 Hours Urine | 4 to 8 °C   | 6 days           |
|                | -20 °C      | 6 months         |
|                | 20 to 25 °C | 2 days           |

**Preparation:**  
**Urine:** work with samples collected in the period of 24 hours. Homogenize the urine, measure its volume and separate 20 mL. Centrifuge for 10 minutes at 3000 rpm and collect the supernatant. Perform a 1:25 dilution of the supernatant with purified water and perform the assay. Multiply the obtained result by 25. If the obtained result surplus the operating range, perform a new dilution by altering the proportion.

### PRODUCT DESCRIPTION

|            |  |   |
|------------|--|---|
| <b>R 1</b> | Picric acid ≥ 10 mmol/L.   |  |
| <b>R 2</b> | Sodium Hydroxide ≥ 200 mmol/L, amino acetic acid ≥ 10 mmol/L, surfactant.  |  |
| <b>STD</b> | Stabilizer; preservative; Creatinine in concentration equivalent to 2 mg/dL. Traceable to the reference material NIST 914. |  |

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

|                                  |           |
|----------------------------------|-----------|
| Autocal H                        | 13.002.00 |
| Normal Control – Quantinorm      | 13.003.00 |
| Pathological Control – Quantialt | 13.004.00 |

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 500 nm (490 – 510 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Incubate the work reagent at 37 °C for 3 minutes.
2. Zero the equipment at 500nm (490 - 510 nm) with purified water.
3. Pipette in the assay tubes:

|              | STD    | Sample |
|--------------|--------|--------|
| STD          | 100 µL | -      |
| Sample       | -      | 100 µL |
| Work Reagent | 1,0 mL | 1,0 mL |

4. Homogenize and insert into the thermostated cuvette at 37 °C. Set the stopwatch.
5. Measure the Sample's and Standard's absorbance at 30 seconds (A1) and 90 seconds (A2).

#### B) CALCULATIONS

Creatinina (mg/dL) =  $\frac{\text{Sample's (A2 - A1)} \times \text{STD Concentration (mg/dL)}}{\text{STD's (A2 - A1)}}$

**Calculations with the Calibration Factor:**  
Calibration Factor =  $\frac{\text{STD Concentration (mg / dL)}}{\text{STD's (A2 - A1)}}$

Creatinina (mg / dL) = Sample's (A2 - A1) x Calibration Factor

**Calculations for urine samples:**  
Creatinina (mg/24 hours) =  $\frac{\text{Creatinina (mg/dL)} \times \text{Urinary Volume (mL)}}{100}$

\*The urine supernatant must be diluted in a 1:25 proportion. To determine the uric acid concentration in the urine sample (mg/dL), multiply the obtained result by 25.

Creatinine (mg/kg/24 hours) =  $\frac{\text{Creatinine (mg/24 hours)}}{\text{Weight (kg)}}$

**Endogenous Creatinine clearance:**  
 Debug (ml/min) =  $\frac{\text{Creatinine in the Urine (mg/dL)} \times \text{Minute Volume}^*}{\text{Creatinine in the serum (mg/dL)}}$

\*Minute Volume = 24 hours urinary volume in mL divided by 1440 minutes.  
 The clearance must be corrected to the patient's body surface.

**Calculation of the Body Surface (SC) without use of the Nomogram:**  
 $SC (m^2) = p^{0,725} \times A^{0,725} \times 0,007184$   
 In which: P = weight in kilograms (kg); A = height in centimeters (cm).

**Clearance:**  
 Clearance (mL/min) =  $\frac{\text{Debug (mL/min)} \times 1,73}{\text{Body Surface (m}^2\text{)}}$

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicalda.ind.br](http://www.biotechnicalda.ind.br)

### C) INTERPRETATION

Creatinine is the product of phosphocreatine decomposition. It is produced endogenously and released into body fluids at a constant rate. Its plasma concentration is maintained within narrow limits, predominantly by glomerular filtration. Consequently, both the plasma creatinine concentration and its renal clearance (creatinine clearance) have been used as markers of glomerular filtration rate.

### PERFORMANCE CHARACTERISTICS

|       | Operating range     |
|-------|---------------------|
| Serum | 0,28 a 12,00 mg/dL  |
| Urine | 3,15 a 321,35 mg/dL |

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

|       | Sensitivity     |                      |
|-------|-----------------|----------------------|
|       | Detection Limit | Quantification Limit |
| Serum | 0,11 mg/dL      | 0,28 mg/dL           |
| Urine | 0,018 mg/dL     | 3,15 mg/dL           |

|       | Analytical Specificity |           |               |
|-------|------------------------|-----------|---------------|
|       | Hemoglobin             | Bilirubin | Triglycerides |
| Serum | 500 mg/dL              | 10 mg/dL  | 500 mg/dL     |
| Urine | 180 mg/dL              | 6 mg/dL   | -             |

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy                    | Serum                | Urine                |
|-----------------------------|----------------------|----------------------|
| Number of Samples           | 40 in duplicate      | 20 in duplicate      |
| Regression Equation         | $y = 0,998x + 0,040$ | $y = 1,002x - 2,343$ |
| Correlation Coefficient (R) | 0,9958               | 0,9995               |

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 1.0 mg/dL and 4.0 mg/dL was 3.8% and 0.8%, respectively.  
**Urine:** by applying the regression equation, the total systematic error estimated for decision levels of 1500 mg/24h and 2000 mg/24h was 0.04% e 0,08%, respectively.

### Precision:

Determined with two runs in duplicate per day for 20 days.

| Serum Samples (mg/dL) | Replicates | Within-Run Precision |        | Total Precision |        |
|-----------------------|------------|----------------------|--------|-----------------|--------|
|                       |            | SD (mg/dL)           | CV (%) | SD (mg/dL)      | CV (%) |
| 0,96                  | 80         | 0,004                | 0,40   | 0,005           | 0,50   |
| 4,07                  | 80         | 0,012                | 0,30   | 0,017           | 0,40   |
| 10,41                 | 80         | 0,125                | 1,20   | 0,132           | 1,30   |

| Urine Samples (mg/dL) | Replicates | Within-Run Precision |        | Total Precision |        |
|-----------------------|------------|----------------------|--------|-----------------|--------|
|                       |            | SD (mg/dL)           | CV (%) | SD (mg/dL)      | CV (%) |
| 93,69                 | 80         | 1,699                | 1,80   | 2,202           | 2,40   |
| 182,10                | 80         | 3,921                | 2,20   | 4,176           | 2,30   |

CV: Coefficient of variation; SD: Standard deviation

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

### REFERENCE RANGES

| Serum / Plasma              | Men       |             | Women     |             |
|-----------------------------|-----------|-------------|-----------|-------------|
|                             | mg/dL     | µmol/L      | mg/dL     | µmol/L      |
|                             | 0,9 - 1,3 | 80 - 115    | 0,6 - 1,1 | 53 - 97     |
| Urine                       | mg/kg/day | µmol/kg/day | mg/kg/day | µmol/kg/day |
|                             | 14 - 26   | 124 - 230   | 11 - 20   | 97 - 177    |
| Clearance (mL/min) - Adults |           | 52,1 - 110  |           |             |

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µmol/L):

Creatinine (mg/dL) x 88,4 = Creatinine (µmol/L)

### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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