Desidrogenase Láctica

Lactate Dehydrogenase (LDH) / Lactato Deshidrogenasa
Ref. 11.004.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine the activity of Lactate Dehydrogenase in serum and plasma samples. Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 30 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

WORKING PRINCIPLE

Method: kinetic UV

Lactate Dehydrogenase (LDH) from sample catalyzes the reduction of pyruvate to lactate in the presence of NADH. The rate of conversion of NADH to NAD⁺ in the reaction is proportional to the catalytic activity of LDH, which can be spectrophotometrically determined by the measurement of the absorbance decrease at 340 nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and heparin plasma.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Serum and Plasma*	4 to 8 °C	1 day
	-20 °C	3 weeks

*Plasma and serum samples must be separated from the cellular elements immediately.

PRODUCT DESCRIPTION

R 1

TRIS buffer ≥ 50 mmol/L; sodium pyruvate ≥ 0,5 mmol/L; activators; stabilizers and preservative.

X

R 2

Carbonate buffer ≥ 2 mmol/L; NADH ≥ 0,5 mM; preservative.

X

Optimized and standardized method according to the recommendations of Deutschen Gesellschaft für Klinische Chemie (DGKC), using calibrated pipettes, manual spectrophotometer and specific absorptivity of the chromogen that provide absolute values.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	REF	13.002.00
Normal Control – Quantinorm		13.003.00
Pathological Control – Quantialt		13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

- Heat the work reagent at 37 °C for 3 minutes.
- Pipette in the test tubes:

	Sample
Sample	20 µL
Work Reagent	1.0 mL

- Homogenize and insert into the thermostated cuvette at 37 °C. Start the stopwatch.
- After 1 minute, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A0 - A1) + (A1 - A2) + (A2 - A3)}{3}$$

The LDH activity is calculate by multiplying the ΔA/min by the factor below:

$$\text{LDH (U/L)} = \Delta A/\text{min} \times 8095$$

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaidta.ind.br

C) INTERPRETATION

LDH corresponds to five isoenzymes, which can be separated according to their electrophoretic mobility. Each isoenzyme is a tetramer composed of two subunits M and/or H. It is widely distributed in the tissues, particularly heart, liver, muscle and kidneys, at concentrations approximately 500 times greater than those normally found in serum. Therefore, it can be used as a marker of acute or chronic tissue injury and, sometimes, in the follow-up of progressive lesions. Serum LDH elevation occurs in several clinical conditions, such as myocardial infarction, hemolysis, neoplasms, liver, kidney, lung and muscle disorders.

PERFORMANCE CHARACTERISTICS

Operating range
11.38 to 2000.00 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosagem and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
5.22 U/L	11.38 U/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
25 mg/dL	40 mg/dL	1500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 0.999x – 6.54
Correlation Coefficient (R)	0.9958

By applying the regression equation, the total systematic error estimated for decision levels of 400 U/L and 1200 U/L was -1.73% and -0.64%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
261.06	80	2.709	1.00	2.834	1.10
776.15	80	2.239	0.30	2.496	0.30
1320.33	80	8.640	0.70	9.191	0.70

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Adults	37 °C
Men	214 - 450 U/L
Women	195 - 453 U/L

Children	37 °C
0 a 6 months	613 - 1020 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

LDH (U/L) x 0,017 = LDH (µKat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
 - The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646
- #### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE
- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
 - To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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