Analytical Specificity			
Bilirubin	Hemoglobin	Triglycerides	Rheumatoid Factor
37.5 mg/dL	500 mg/dL	500 mg/dL	60 IU/mL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y = 1.0099 + 0.0038
Correlation Coefficient (R)	0.9936

Precision:
Within-run precision was determined with 3 lots and 10 replicates in an analytical run; between-run precision was determined with 3 lots and 3 replicates in an analytical run.

Samples (mg/L)	Within-Run Precision		Between-Run Precision	
	SD (mg/L)	CV (%)	SD (mg/L)	CV (%)
2.203	0.083	3.78	0.073	2.76

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices’ instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Serum and Plasma	≤ 0.5 mg/L
------------------	------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Unit conversion:
D-dimer (mg/L) x 1000 = D-dimer (ng/mL)

D-dimer results can also be presented in FEU (Fibrinogen Equivalent Units):
D-dimer (ng/mL) x 1,74 = D-dimer (ng/mL FEU).

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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