Electrolyte IFS Solution

Electrolyte IFS Solution/Electrolyte IFS
Solución
Ref. 40.021.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Product used to complement selective ion electrodes that perform quantitative determinations of sodium, potassium, chloride and ionized calcium. Use for *in vitro* diagnostics only.

STORAGE AND HANDLING CONDITIONS

- Store at 15 to 30 °C. Keep away from light.
- Expiry date stated on the label.
- Do not use products whose expiration date has passed.

WORKING PRINCIPLE

Product used to complement selective ion electrodes that perform quantitative determinations of sodium, potassium, chloride and ionized calcium.

PRODUCT DESCRIPTION

NaHEPES 5%, HEPES 5%, sodium chloride 10%, potassium chloride 5%, calcium chloride 2%, lithium acetate 0.01%, preservative.

INSTRUCTION FOR USE

- Shake gently before use.
- Fill 3/4 of the electrodes with IFS Solution.
- Allow electrodes to stabilize for 5-10 minutes.

RESIDUAL RISKS, CARE AND PRECAUTIONS

- Use PPE and carry out procedures in accordance with Good Clinical Laboratory Practices.
- Follow the requirements recommended in Good Clinical Laboratory Practices for water used in the Laboratory.
- Do not use the product when it presents visual characteristics that do not comply with those specified in the MSDS.

WARNINGS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information is described in the Individual Chemical Product Safety Data Sheet (MSDS) for this product, available at www.biotecnica.ind.br or by calling +55 35 3214-4646.
- Dispose of reaction residues in accordance with Good Clinical Laboratory Practices and the Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / SAC - CUSTOMER CARE SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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