# Electrolyte ISE Pak

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**ENGLISH**

**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**
**INTENDED USE**

 Product intended for calibration of analytical tests in Ion-Selective Electrode (ISE) methodology for quantitative determination of sodium, potassium, chloride and ionized calcium. Use for *in vitro* diagnostics only.

**STORAGE AND HANDLING CONDITIONS**

- Store at 15 to 30 °C. Keep away from light.
- Expiry date stated on the label.
- Once opened, the product in use is stable for 7 months, as long as the recommended storage conditions are followed.
- Do not use products whose expiration date has passed.

**WORKING PRINCIPLE**

Ion-selective electrodes are devices used in laboratory analyzes to measure the amount of ions in biological samples, such as blood, serum, plasma and urine. Each ion selective electrode is designed to detect a single type of ion, such as sodium, potassium, chloride or calcium. The operating principle of these electrodes is based on the selective ability of the sensor to respond only to a specific type of ion. Each electrode has a selective membrane that allows only specific ions to pass through it. When these ions cross the membrane, an ion exchange occurs that results in the generation of a voltage or membrane potential. This potential is measured between the sample and the electrode membrane and the ion concentration can be determined using the Nernst equation.

**PRODUCT DESCRIPTION**
**SOLUTION A**

 Sodium: 150,0 mmol/L  
 Potassium: 5,0 mmol/L  
 Chloride: 115,0 mmol/L  
 Ionized calcium: 0,9 mmol/L

**SOLUTION B**

 Sodium dio: 100,0 mmol/L  
 Potassium: 1,8 mmol/L  
 Chloride: 72,0 mmol/L  
 Ionized calcium: 1,5 mmol/L

**REFERENCE SOLUTION**

 Potassium Chloride: 8,9%  
 Preservative > 0,05%

**WASTE CONTAINER**

**INSTRUCTION FOR USE**

- Shake gently before use.
- Before use, remove the red plastic cover from the connector. Keep this cap to close the connectors on used packaging.
- Write the installation date on the label side of the ISE Pak.
- Connect the ISE Pak to the left side of the analyzer.
- The use of this solution must be in accordance with what is specified in the device manufacturer's manual.

**PERFORMANCE CHARACTERISTICS**

Operational range:

Parameters	Whole blood/Serum/Plasma/Aqueous solutions
Sodium	40 – 205 mmol/L
Potassium	1,5 – 15,0 mmol/L
Chloride	50 – 200 mmol/L
Ionized calcium	0,2 – 5,0 mmol/L

**REFERENCE INTERVALS**

Parameters	Whole blood/Serum/Plasma
------------	--------------------------

Sodium	136 – 145 mmol/L
Potassium	3,5 – 5,1 mmol/L
Chloride	97 – 108 mmol/L
Ionized calcium	1,15 – 1,35 mmol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

### RESIDUAL RISKS, CARE AND PRECAUTIONS

- Contains toxic preservatives that have potential to irritate the skin. Therefore, if the reagent comes into contact with skin or clothing, rinse immediately with plenty of water and consult a doctor if skin irritation develops.
- Use PPE and carry out procedures in accordance with Good Clinical Laboratory Practices.
  - Follow the requirements recommended in Good Clinical Laboratory Practices for water used in the Laboratory.
  - Do not use the product when it presents a visual characteristic that does not comply with what is specified in its MSDS.

### WARNINGS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information is described in the Individual Chemical Product Safety Data Sheet (MSDS) for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by calling +55 35 3214-4646.
- Dispose of reaction residues in accordance with Good Clinical Laboratory Practices and the Health Service Waste Management Program (PGRSS).

### QUALITY ASSURANCE / SAC - CUSTOMER CARE SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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