FR Látex
 RF Latex / FR Látex
 Ref. 30.004.00 (KIT) / 30.004.00 (REAGENTE)
 Responsável Técnico:
 Dr. Gilson Serio Pizzo
 CRF MG – 5310
 Anvisa 80027310270

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the determination Rheumatoid Factor (RF) in serum. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagents and controls.

WORKING PRINCIPLE

Method: Plate agglutination

RF Latex is a plate agglutination method for qualitative and semi-quantitative detection, where latex particles coated with human gamma globulin are agglutinated by rheumatoid factors present in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	8 days
	-20 °C	3 months

PRODUCT DESCRIPTION

Suspension of Latex particles coated with human gamma globulin, pH 8.2. Preservative. The sensitivity of RLAX is traceable to NIBSC 64/002 (National Institute for Biological Standards and Control).

CONTROL + Human Serum with FR concentration > 30 UI/mL. Preservative. (INCLUDED IN THE RF LATEX KIT)

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CONTROL -

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

CONTROL (+) (INCLUDED IN THE RF LATEX KIT)	REF	30.004.00
CONTROL (-) (INCLUDED IN THE RF LATEX KIT)		30.004.00

NECESSARY EQUIPMENT FOR TESTING

- Reaction plate and disposable rod;
- Pipette, clock or stopwatch..

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Automação: This procedure cannot be applied to automated analyzers.

Qualitative Method

- Before starting the tests, allow the RLAX, controls and samples to reach room temperature. Homogenize before use.
- On the reaction plate use a different circle for each sample or control.
- Place 20 µL of the sample / controls.
- In each circle add 20 µL of RLAX.
- Mix each determination using disposable rods (one for each sample).
- Homogenize the plate for 2 minutes in circular motions, tilting it back and forth. Observe the result under artificial light: Presence or absence of agglutination.

Semi-quantitative method

- Dilute the sample with saline (0.9% NaCl) as follows:

Dilution	Saline (µL)	Sample(µL)	Dilution Factor
1/2	100	100 (pure)	2
1/4	100	100 (sol.1/2)	4
1/8	100	100 (sol.1/4)	8
1/16	100	100 (sol.1/8)	16

- Carry out the test of each dilution according to the qualitative method.

B) CALCULATIONS

The approximate concentration in the semi-quantitative test is obtained by multiplying the dilution factor by the sensitivity of the test (8 UI/mL).

Example:

Agglutination up to 1/8 dilution.

Dilution Factor = 8

FR = 64 UI / mL (8 x 8)

C) INTERPRETATION
 Rheumatoid factors are a group of antibodies that bind to the Fc fraction of immunoglobulins G. Although they are present in a large number of rheumatic disorders, such as systemic lupus erythematosus and Sjögren's syndrome, their main clinical interest is in the diagnosis of rheumatoid arthritis..

PERFORMANCE CHARACTERISTICS

Operating range
Up to 600 UI/mL

Sensitivity
Detection Limit
8 UI/mL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
1000 mg/dL	40 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	100
Positive Agreement	100%
Negative Agreement	96%

Prozone High Dose Effect:

There is no prozone effect up to 600 IU/mL of FR.

Precision:

Intra-assay: the study was carried out in two daily runs, in duplicate, for 10 days, using 04 samples, two positive and two negatives, with 100% agreement.

Inter-assay: the study was carried out in a daily run of four replicates, for 5 days, using 04 samples, two positive and two negatives, with 100% agreement

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Adults	Up to 8 UI/mL
--------	---------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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• SHMERLING, Robert H.; DELBANCO, Thomas L. The Reumatoid Factor: An Analysis of Clinical Utility. *The American Journal Of Medicine*, v. 91, p.528-534, 1991.

• WOLFE, Frederick; CATHEY, Mary Ann; ROBERTS, F. Kay. The Latex Test Revisited: Rheumatoid Factor Testing in 8,287 Rheumatic Disease Patients. *Arthritis And Rheumatism*, v. 34, n° 8, p.951-960, 1991.

• WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

• YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AAC Press, 2000

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