Fator Reumatóide Turbilátex

Rheumatoid Factor Turbilatex / Factor Reumatoide Turbilátex

Ref. 20.007.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for the determination of Rheumatoid Factors (RF) in serum. *In vitro* diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- **RDIL e RLAX:** After being opened, reagents in use are stable up to the expiration date printed on the label, provided that they are followed at recommended storage conditions (2 to 8 ° C). Ready-to-use reagents.
- Do not use reagents whose shelf life has expired.
- Lyophilized Calibrator: Reconstitute with 2.0 mL of purified water. Leave to stand for 10 minutes. Gently homogenize by inversion avoiding foaming. After reconstitution is stable for 30 days if stored at 2 to 8 ° C

WORKING PRINCIPLE

Method: Immunoturbidimetric

latex particles coated with human gamma globulin are agglutinated by RF present in the sample. The agglutination causes a turbidity whose intensity is proportional to the concentration of RF.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

	Temperature	Stability Period
	2 to 8 °C -20 °C	7 days 3 months

PRODUCT DESCRIPTION

R

DIL

Tris buffer 20 mmol / L pH 8.2, preservative.

R

LAX

Latex particles coated with gammaglobulin pH 7.4; preservative.

CAL

Human serum containing Rheumatoid Factor; preservative.
The RF value is traceable to the Rheumatoid Arthritis Serum NIBSC Code: 64/2 WHO.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Rheumatoid Factor Turbilatex CAL (Kit Included)	REF	20.007.00
Rheumatic Control Level 1		21.007.00
Rheumatic Control Level 2		21.008.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 650 nm (600 - 550 nm).
- Glass and / or automatic pipettes.
- Test tubes.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaInd.br

REACTION PARAMETERS	
Type of reaction	Fixed time
Wave length	650 nm (600 – 650 nm)
Temperature	37 °C
Sample/CAL*	7 µL
RDIL*	800 µL
RLAX*	200 µL
Reading 1 (A1)	Immediately
Reading 2 (A2)	2 minutes
Calibration	Calibration curve**

Sample volumes/CAL, RDIL and RLAX be modified as long as the established ratio is maintained.

****Calibration Curve:**

Dilute the CAL in NaCl 9g / L according to follow:

Dilution	1	2	3	4	5	6
CAL (µL)	-	25	50	100	200	400
NaCl 9g/L (µL)	400	375	350	300	200	-
Factor	0	0,0625	0,125	0,25	0,5	1

Multiply the concentration of Rheumatoid Factor Turbilatex Turbilátex CAL, indicated on the label, by the corresponding factor to obtain the RF concentration of each dilution.

B) CALCULATIONS

The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION

Rheumatoid factors are a group of antibodies that bind to the Fc fraction of immunoglobulins G. Although they are present in a large number of rheumatic disorders, such as systemic lupus erythematosus and Sjögren's syndrome, their main clinical interest is in the diagnosis of rheumatoid arthritis.

PERFORMANCE CHARACTERISTICS

Operating range
10.3 to 141 UI/mL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
3.7 UI/mL	10.3 UI/mL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
100 mg/dL	18 mg/dL	750 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	38 in duplicate
Regression Equation	y = 0.993x - 0.197
Correlation Coefficient (R)	0.9986

By applying the regression equation, the total systematic error estimated for decision levels of 50 UI/mL and 100 UI/mL was -1.1% and -0.9%, respectively.

Prozone High Dose Effect:

There is no prozone effect up to 650 UI/mL of FR.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (UI/mL)	Replicates	Within-Run Precision		Total Precision	
		SD (UI/mL)	CV (%)	SD (UI/mL)	CV (%)
18,4	80	0,344	1,9	0,401	2,2
47,2	80	1,152	2,4	1,449	3,1
112,8	80	1,414	1,3	2,133	1,9

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS
- Use protective equipment in accordance with the Good Laboratory Practices. - Follow the Good Laboratory Practices' instructions to establish the quality of water. - Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS. - Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES	
Adults	Up to 20 UI/mL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material. - The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE
- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br. - To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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