Ferritina

Ferritin | Ferritina

Ref. 20.008.00

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Dr. Gilson Serio Pizzo

CRF MG - 5310

Anvisa 80027310249

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ENGLISH
 **BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

- INTENDED USE**
Kit intended to determine ferritin in serum samples. Diagnostic use only.
- STORAGE AND HANDLING**
- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
 - Reagents ready for use. Lyophilized calibrator.
 - Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
 - Do not use reagents whose shelf life has expired.

PREPARATION AND HANDLING
Calibrator: reconstitute with 3,0 mL of purified water. Wait 10 minutes at room temperature. Homogenize gently by inversion, avoiding foam formation. Reconstituted, the calibrator is stable for 30 days when stored at 2 to 8 °C, and 3 months when stored at -20 °C.

WORKING PRINCIPLE
Method: Immunoturbidimetric.
Latex particles coated with anti-human ferritin antibodies form insoluble complexes when mixed with the ferritin present in the sample. The turbidity intensity is proportional to the ferritin concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND PRESERVATION
Sample Type: serum.
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

Serum	Temperature	Sample Stability
	2 to 8 °C -20 °C	7 days 3 months

R

DIL

R

LAX

CAL

TRIS Buffer > 20 mmol/L pH 8.2; preservative.

Latex particles coated with anti-Ferritin antibody pH 8.2; preservative

Human serum containing Ferritin. The value of Ferritin is traceable to the 2nd International Standard for Ferritin (80/578) 1992 (WHO) NIBSC.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and control below:

Ferritin CAL (included in the kit)

Medium Multi-parameter Control

REF

20.008.00

21.003.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 540 nm.
- Glass and / or automatic pipettes and test tubes.

TEST PROCEDURE, CALCULATION, AND INTERPRETATION
A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br.

Reaction Type	Fixed Time
Wavelength	540 nm
Temperature	37 °C
Sample Volume / CAL*	18 µL
RDIL Volume*	160 µL
RLAX Volume*	40 µL
First Reading	Immediately
Second Reading	5 minutes
Calibration**	Calibration Curve

* Volumes can be modified as long as the established proportion is maintained.

****CALIBRATION CURVE:**

Dilute the Ferritin CAL in a 9,0 g/L NaCl solution as follows:

Dilution	1	2	3	4	5	6
Calibration (µL)	-	15	25	50	100	200
NaCl 9,0 g/L (µL)	200	225	175	150	100	-
Fator	0	0,0625	0,125	0,25	0,5	1

Multiply the concentration of the Ferritin CAL indicated on the label by the corresponding dilution factor to obtain the Ferritin concentration of each dilution.

B) CALCULATIONS
The equipment automatically calculates the analyte concentration in each sample.

C) INTERPRETATION
Ferritin is found in almost every cell in the body. In hepatocytes and macrophages, ferritin provides an iron reserve readily available for the formation of hemoglobin and other hemeproteins. It is stored so that the iron is protected from body fluids and is unable to cause oxidative damage, which could happen if it were in its ionic form. Small amounts of ferritin are present in serum at concentrations approximately proportional to the total body iron storage. On the other hand, liver damage and many pathological processes result in the release of relatively large amounts of ferritin in the plasma. Ferritin can be considered the most sensitive marker for detecting the early stages of iron deficiency when compared to hemoglobin levels, especially after 6 months of age, when children are more vulnerable to iron deficiency anemia due to depletion of reserves from pregnancy and low dietary iron intake.

PERFORMANCE CHARACTERISTICS

Operating range		
9,9 µg/L to 585,0 µg/L		

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
2.7 µg/L	9.9 µg/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
500 mg/dL	35 mg/dL	250 mg/dL

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 1.002x + 0.405
Correlation Coefficient (R)	0.9999

By applying the regression equation, the total systematic error estimated for decision levels of 100 µg/L and 200 µg/L was 0.6% and 0.4%, respectively.

Prozone High Dose Effect:
Prozone effect does not occur up to 2300 µg/L of ferritin.

Precision: Determined with two runs in duplicate per day for 20 days.

Samples (µg/L)	Replicates	Within-Run Precision		Total Precision	
		SD (µg/L)	CV (%)	SD (µg/L)	CV (%)
154.80	80	3.19	2.10	4.52	2.90
264.90	80	4.84	1.80	5.76	2.20
407.30	80	5.80	1.40	5.87	1.40

% CV: Coefficient of variation; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices’ instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Children and Teenagers	
Newborns	25 – 200 µg/L
1 month	200 – 600 µg/L
2 - 5 months	50 – 200 µg/L
6 months - 15 years	7 – 140 µg/L

Adults	
Men	20 – 250 µg/L
Women	10 – 120 µg/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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