Ferro Ferrozine

Iron Ferrozine / Hierro Ferrozina
Ref. 12.005.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine iron in serum samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (15 to 30 °C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Ferrozine

The iron from the sample, at acidic pH, is released from transferrin. Fe³⁺ ions are reduced by hydroxylamine to Fe²⁺. These react with the chromogen ferrozine to form a colored complex, which can be spectrophotometrically determined at 560 nm. Color intensity is proportional to the concentration of iron in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Sample Stability
Serum	4 to 8 °C	3 weeks
	20 to 25 °C	7 days
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1	Succinate buffer > 0.1 mol/L; hydroxylamine hydrochloride > 0.1 mol/L; detergent; reducing agent.	X
R 2	Ferrozine > 1.0 mmol /L; preservative.	
STD	Solution of Fe ³⁺ ions in concentration equivalent to 100 µg/dL, preservative, stabilizer. Traceable to reference materials NIST-3126a and NIST-928.	

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Autocal H	13.002.00	REF
Normal Control – Quantinorm	13.003.00	
Pathological Control – Quantial	13.004.00	

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 560 nm (540 – 580 nm).

- Test tubes, glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Reagent's Blank (B1)	Standard (P1)	Sample's Blank (A1)	Sample (A2)
Purified water	450 µL	250 µL	-	200 µL	-
Sample	-	-	-	250 µL	250 µL
STD	-	-	250 µL	-	-
R1	800 µL	800 µL	800 µL	800 µL	800 µL
R2	-	200 µL	200 µL	-	200 µL

2. Homogenize and incubate at 15 to 30 °C for 10 minutes.
3. Measure the Reagent's Blank (B1), Sample's Blank (A1), Sample's (A2) and Standard's (P1) absorbance, setting the equipment to zero with the Blank. Measure the absorbance at 560 nm (540 – 580 nm).

B) CALCULATIONS

Iron (µg/dL) = $\frac{A2 - A1 - B1}{P1} \times \text{STD Concentration (µg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (µg/dL)}}{(P1 - B1)}$

Iron (µg/dL) = (A2 - A1 - B1) x Calibration Factor

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicaltda.ind.br

C) INTERPRETATION

The absorption of iron ingested as Fe²⁺ happens in the duodenum and upper jejunum. The free trivalent form and Fe³⁺ bound to the hemoproteins, present in the diet, are reduced by vitamin C. About 1 mg of iron is absorbed daily. Once in the mucosal cells, Fe²⁺ ions bind to transporters. Before passing to the plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bind to apotransferrin. The apotransferrin-Fe³⁺ complex is called transferrin, which carries the iron ions in the blood plasma. The determination of iron (non-heme) is used in the diagnosis and treatment of diseases such as anemia due to iron deficiency, hemochromatosis, chronic kidney disease and in the monitoring of other types of anemias.

PERFORMANCE CHARACTERISTICS

Operating range
6.66 µg/dL to 500.0 µg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
4.24 µg/dL	6.66 µg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	40 mg/dL	900 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1.002x - 0.165$
Correlation Coefficient (R)	0.9999

By applying the regression equation, the total systematic error estimated for decision levels of 100 µg/dL and 350 µg/dL was 0.03% and 0.15%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (µg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (µg/dL)	CV (%)	SD (µg/dL)	CV (%)
64.78	80	0.711	1.10	1.932	3.00
222.30	80	3.145	1.40	6.092	2.70
418.84	80	1.233	0.30	6.640	1.60

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory shall establish the chemical, microbiological and particle requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been

- defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.
 - Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
 - Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Men	65 – 170 µg/dL
Women	50 – 170 µg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µmol/L):

Iron (µg/dL) x 0,179 = Iron (µmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotecnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnicos products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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- STOOKEY, L. L. Ferrozina - a new spectrophotometric reagent for iron. Anal. Chem. v.42, p.779-781, 1970.
- ARTISS, J.D.; VINOGRADOV, S.; ZAK, B. Spectrophotometric study of several sensitive reagents for serum iron. Clin. Biochem. v.14, p.311-315, 1981
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.
- BURTIS, Carl A.; ASHWOOD, Edward R. Tietz: Fundamentos de Química Clínica. 4. ed. Philadelphia: W. B. Saunders, 1996. 836 p.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES			
	Consultar as instruções para utilização Consult instructions for use Consúltense las instrucciones de uso		Descartar corretamente Dispose properly Desechar adecuadamente
REF	Número de catálogo Catalog number Número de catálogo	R <N>	Reagente Reagent Reactivo
LOT	Código do lote Batch code Código de lote		Límite de temperatura Temperature limitation Límite de temperatura