PERFORMANCE CHARACTERISTICS

Operating Range:

The operating range of the product is delimited by the concentration of the 1/5 and 1/40 dilutions of the calibration curve. For values outside the operating range, dilute the sample with the Fibrinogen RDIL (buffer), proceed with a new dosage and multiply the result by the dilution factor. Example

- If the clotting time of the sample is greater than the clotting time of diluting 1/40 of the calibration curve, repeat the test diluting the sample in the proportion of 1/5. Multiply the obtained result by 5 to determine the concentration of fibrinogen in the sample. Do not test dilutions lower than 1/3.
- If the clotting time of the sample is lower than the clotting time of the 1/5 dilution of the calibration curve, repeat the test diluting the sample in the proportion of 1/20. Multiply the obtained result by 20 to determine the concentration of fibrinogen in the sample.

Sensitivity:

The sensitivity of the product is defined by the concentration of the 1/40 dilution of the calibration curve.

Analytical Specificity

Hemoglobin

400 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy

Number of Samples

50

Regression Equation

y = 1.1453 - 23,751

Correlation Coefficient

0.9944

Precision:

Within-run precision was determined with 3 lots and 5 replicates for 5 days; between-run precision was determined with 3 lots and 5 replicates for 5 days.

Samples (mg/dL)

Within-Run Precision

Between-Run Precision

SD (mg/dL)

CV (%)

SD (mg/dL)

CV (%)

136.60

6.77

4.93

1.70

1.30

274.20

4.50

1.67

1.30

0.50

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Plasma

200 to 400 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges. To convert to g/L, multiply by the factor 0.01

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. **Washington DC: AACC Press**, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.
- Clauss A. Acta Haematol 1957;17:237–246
- International commitee For Standardization In **Haematology And International Committee On Trombosis And Haemostasis**. Amer. J. Clin. Path. :88: 779 1985.
- WORLD HEALTH ORGANIZATION et al. **Use of anticoagulants in diagnostic laboratory investigations**. World Health Organization, p. 34, 2002.
- NCCLS. **Procedure for the Determination of Fibrinogen in Plasma; Approved Guideline—Second Edition**. Clinical and Laboratory Standards Institute, v. 21, n. 18, 2001.
- NCCLS. **Procedure for the Determination of Fibrinogen in Plasma; Approved Guideline—Second Edition**. Clinical and Laboratory Standards Institute, v. 21, n. 18, 2001.
- Biotécnica: Dados de arquivo

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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In Vitro Diagnostic medical device

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