## FLUORETO

FLUORIDE / FLUORURO  
Ref. 35.003.00

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### ENGLISH

 BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

### INTENDED USE

Anticoagulant solution to obtain plasma for glucose determination. Use *in vitro* diagnostic.

### STORAGE AND HANDLING

- Store at 15 to 30 °C and protect from light.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagent.

### WORKING PRINCIPLE

EDTA blocks the ionized calcium by chelation. Thus, the coagulation process (Calcium dependent) is paralyzed, inhibiting the conversion of Prothrombin to Thrombin and the action of Thrombin on Fibrinogen that forms Fibrin (clot). Fluoride inhibits glycolysis.

### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

### PRODUCT DESCRIPTION

**R 1** Potassium fluoride > 1 mol/L; EDTA 100 mmol/L and light green > 10 µmol/L. 

### NECESSARY EQUIPMENT FOR TESTING

- Test tubes.
- Centrifuge.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Use 01 drop of anticoagulant for every 3 mL of blood.
2. Homogenize gently by inversion.
3. Centrifuge at 3000 rpm for 5 minutes and separate the plasma.

Note: Anticoagulant / blood volumes may be modified provided that the established ratio is maintained

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.

- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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