2.45 U/L

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	300 mg/dL	40 mg/dL	750 mg/dL
Felinos	500 mg/dL	40 mg/dL	750 mg/dL
Equinos	200 mg/dL	30 mg/dL	2000 mg/dL
Bovinos	500 mg/dL	40 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	40 in duplicate
Regression Equation	y = 0.998x - 0.255
Correlation Coefficient (R)	1.0000
Accuracy - Felinos	
Number of Samples	20 in duplicate
Regression Equation	y = 0.990x + 0.705
Correlation Coefficient (R)	0.9999
Accuracy - Equinos	
Number of Samples	11 in duplicate
Regression Equation	y = 1.002x - 3.087
Correlation Coefficient (R)	0.9972
Accuracy - Bovinos	
Number of Samples	10 in duplicate
Regression Equation	y = 1.032x - 0.4286
Correlation Coefficient (R)	0.9993

Precision:

Determined with two runs and in 20 replicates

Samples (U/L)	Number of samples	Within-Run Precision		Samples (U/L)	Number of samples	Between-Run Precision	
		SD (U/L)	CV (%)			SD (U/L)	CV (%)
146.24	20	6.79	4.60	143.40	40	6.70	4.70
482.19	20	12.34	2.60	482.48	40	11.98	2.50

%CV: Coefficient of Variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Species	37 °C
Canine	20.0 - 156.0 U/L
Feline	25.0 - 93.0 U/L
Equine	143.0 - 395.0 U/L
Bovine	0.0 - 488.0 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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INTERVALO DE REFERENCIA

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Caninos	20,0 - 156,0 U/L
Felinos	25,0 - 93,0 U/L
Equinos	143,0 - 395,0 U/L
Bovinos	0,0 - 488,0 U/L

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- YOUNG, D.S. *Effects of drugs on clinical laboratory tests* - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- KANEKO, Jiro J; HARVEY, John W; BRUSS, Michael L. *Clinical biochemistry of domestic animals*. 5th. ed. San Diego: Academic Press, c1997 932p.
- THRALL, Mary Anna. *Hematologia e bioquímica clínica veterinária*. São Paulo: Roca, 2007. 582 p.

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