Fosfatase Alcalina - IFCC

Alkaline Phosphatase – IFCC /
Fosfatasa Alcalina - IFCC
Ref. 11.065.00

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ENGLISH


⚠ BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE
Kit for the quantitative determination of the enzymatic activity of Alkaline Phosphatase in serum and plasma. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.
- Reagents ready for use. The reagents can also be used as a Work Reagent.
- Work Reagent:** mix in the proportion of 1 part of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 14 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).

WORKING PRINCIPLE
Method: Kinetic / adapted IFCC
The determination is based on a kinetic and colorimetric method, in which the sample alkaline phosphatase, in an alkaline medium, hydrolyzes the p-nitrophenylphosphate present in the reagent, releasing p-nitrophenol and inorganic phosphate. The concentration of p-nitrophenol, measured at 405 nm, is directly proportional to the enzymatic activity of alkaline phosphatase in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY
Sample Type: serum and heparin plasma.
Collection, handling and preparation: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.
Preservation:

	Temperature	Stability Period
Serum and Plasma	4 a 8 °C -20 °C	7 days 2 months

PRODUCT DESCRIPTION

R 1	Buffer, ativadores, 2-Amino 2-Methyl-1-Propanol < 800 mmol/L	
R 2	p-Nitrophenylphosphate < 50 mmol/L, preservative	X

QUALITY CONTROL
The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it's recommended the use of the calibrator and controls below:

Serum Calibrator - Autocal H		13.002.00
Normal Control Serum - Quantinorm	REF	13.003.00
Pathological Control Serum - Quantial		13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm (390 – 425 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE
1. Pipette in the test tubes:

	Bi-reactant protocol			Monoreagent Protocol		
	Blank	Calibrator	Sample	Blank	Calibrator	Sample
R1	500 µL	500 µL	500 µL	RT	1000 µL	1000 µL
Incubate at 37°C for 180 seconds						
R2	500 µL	500 µL	500 µL	CAL	-	25 µL
CAL	-	25 µL	-	Sample	-	25 µL
Sample	-	-	25 µL			

- Homogenize and incubate the reaction tubes at 37 °C for 54 seconds;
- Immediately after incubation, measure the initial absorbance (A₀) of each reaction at 405 nm;
- Take new readings every 60 seconds (A₁ and A₂).

B) CALCULATIONS
Alkaline phosphatase (U/L) = Calibration Factor * [(ΔA_{sample}/ min) – (ΔA_{blank}/min)]

In what:

$$\Delta A / \text{min} = \frac{(A_1 - A_0) + (A_2 - A_1)}{2}$$

Calibration Factor = $\frac{\text{Calibrator Concentration (U/L)}}{\Delta A_{\text{calibrador}} / \text{min} - \Delta A_{\text{blank}} / \text{min}}$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicalda.ind.br

C) INTERPRETATION
The enzyme Alkaline phosphatase is present in the body in membranes and cell surfaces located in the mucosa of the small intestine and in the proximal contorted tubules of the kidney, bone (osteoblasts), liver and placenta. Although its exact metabolic function is not fully understood, it is speculated that alkaline phosphatase is associated with the transport of lipids in the intestine and with the process of bone calcification. The enzyme is elevated in different clinical conditions: cholestasis, diseases of the skeletal system (such as Paget's disease), hyperparathyroidism, rickets, osteomalacia, as well as cases of fractures and malignant tumors. The physiological increase may occur in children and young people as a result of increased osteoblastic activity due bone development.

PERFORMANCE CHARACTERISTICS

Operating range
6.5 to 1600 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0,55 U/L	4,09 U/L

Analytical Specificity			
Hemoglobin	Bilirubin	Triglycerides	Ascorbic Acid
250 mg/dL	30 mg/dL	2000 mg/dL	100 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of samples	40 in duplicate
Regression equation	y = 0.9511x – 0.6514
Correlation coefficient (R)	0,9999

By applying the regression equation, the total systematic error estimated for decision levels of 50 U/L, 150 U/L e 400 U/L was -6,2%, -5,3% e -5,1%, respectively.

Precision:
Within-run precision was determined with 5 replicates for 10 days; between-run precision was determined in 3 conditions, with 3 replicates for 5 days.

Samples (U/L)	Within-Run Precision		Samples (U/L)	Between-Run Precision	
	SD (mg/dL)	CV (%)		SD (U/L)	CV (%)
72.10	1.55	2.1%	71.28	1.35	1.9%
188.34	3.45	1.8%	190.93	5.13	2.7%
330.31	1.98	0.6%	332.33	3.21	1.0%

%CV: Coefficient of Variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

0 - 6 years old	81,9 – 350,3 U/L
6 - 12 years old	48,8 – 445,9 U/L
Male	40 – 130 U/L
Female	35 – 105 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotecnica.ind.br.

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- BURTIS, C. A.; ASHWOOD, E. R.; BRUNS, D. E. **Tietz Fundamentos de Química Clínica**, Saunders Elsevier, 6 ed, 2008.
- Gu Wanjiang, MD, Hu Jie, MB, Ge Liang, MB, Wang Cheng, MD, Xu Tian, MB, Shen Jianjiang, MB, Zhang Chunni, MD, Establishment of Reference Interval for Alkaline Phosphatase in Healthy Children of Various Ethnicities, Aged 0-12 Years, Laboratory Medicine, Volume 48, Issue 2, May 2017, Pages 166–171.
- Thomas L, Müller M, Schumann G, Weidemann G et al. Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. J Lab Med 2005;29:301-308.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES

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