Fosfatase Alcalina - DGKC

Alkaline Phosphatase-DGKC / Fosfatasa Alcalina-DGKC
Ref. 11.064.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine the activity of Alkaline Phosphatase in serum and plasma samples. Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- **Work Reagent:** mix in the proportion of 4 parts of R1 + 1 part of R2. Homogenize gently. The work reagent is stable for 28 days, as long as the preparation instructions and the recommended storage conditions are followed (2 to 8 °C).
- The reagents can be used in a bi-reactant mode, as long as the established proportions are followed. Consult the Scientific Advisory.

WORKING PRINCIPLE

Method: DGKC kinetics

Alkaline Phosphatase from sample, at alkaline pH, hydrolyzes p-nitrophenylphosphate, releasing p-nitrophenol and phosphate. The catalytic concentration of the enzyme present in the sample is determined from the rate of formation of p-Nitrophenol, measured at 405nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and heparin plasma

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	7 days
	-20 °C	2 months

PRODUCT DESCRIPTION

R 1 Diethanolamine ≥ 5,0% v/v; magnesium chloride ≥ 0,5 mmol/L; preservative. 

R 2 Disodium PNPP ≥ 20 mmol/L; stabilizer; preservative. 

Traceable to the specific absorbance of the chromogen p-nitrophenol using calibrated pipettes and handheld spectrophotometer that provide absolute values, based on modified methodology of the DGKC (German Society of Clinical Chemistry).

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantial	13.004.00

REF

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Heat the work reagent at 37 °C for 3 minutes.
2. Zero the equipment at 405nm with purified water.
3. Pipette in the test tubes:

	Volume
Sample	20 µL
Work Reagent	1,0 mL

4. Homogenize and insert in the thermostated cuvette at 37 °C. Activate the chronometer.
5. After 1 minute, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

The Alkaline Phosphatase activity is calculate by multiplying the ΔA/min by the factor below:

$$\text{Alkaline Phosphatase (U/L)} = \Delta A/\text{min} \times 2720$$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicalda.ind.br

C) INTERPRETATION

Alkaline phosphatase activity is present in the body in membranes and cell surfaces located in the mucosa of the small intestine and in the proximal contorted tubules of the kidney, bone (osteoblasts), liver and placenta. Although its exact metabolic function is not understood, it appears that alkaline phosphatase is associated with the transport of lipids in the intestine and with the process of bone calcification. The enzyme is elevated in different clinical conditions: cholelstatias, diseases of the skeletal system (such as Paget's disease), hyperparathyroidism, rickets, osteomalacia, as well as cases of fractures and malignant tumors. Physiological increase may occur in children and young people as a result of increased osteoblastic activity in accelerated bone growth.

PERFORMANCE CHARACTERISTICS

Operating range
2,41 a 700,00 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1,68 U/L	2,41 U/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
200 mg/dL	40 mg/dL	500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 1.001x – 1.20
Correlation Coefficient (R)	0.9996

By applying the regression equation, the total systematic error estimated for decision levels of 200 U/L and 500 U/L was -0.50% and -0.14%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
174.10	80	0.865	0.50	0.916	0.50
371.61	80	0.626	0.20	1.329	0.40
633.23	80	0.704	0.10	0.704	0.10

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.

- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	37 °C
Men	< 270 U/L
Women	< 240 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µkat/L):

Alkaline Phosphatase (U/L) x 0,017 = Alkaline Phosphatase (µkat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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