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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.


INTENDED USE
Kit intended to determine phosphorus acid in serum and urine samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed).
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Daly and Ertingshausen modified.
Inorganic phosphate from sample reacts with ammonium molybdate in an acidic medium, producing the phosphomolybdate complex. The latter is spectrophotometrically measured at 340 nm. Color intensity is proportional to the inorganic phosphate concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: serum and urine.
Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum	4 to 8 °C	4 days
	20 to 25 °C	1 day
	-20 °C	1 year
24 Hours Urine	4 to 8 °C	6 months (in pH < 5,0)
	20 to 25 °C	2 days

Preparation:
Urine: work with samples collected in the period of 24 hours in a flask containing 20 to 30 mL of HCl 50% to avoid phosphate precipitation. Homogenize the urine, measure its volume and separate 20 mL. Adjust the pH below 3,0. In case of precipitation, centrifuge the sample and utilize the supernatant. Perform a 1:10 dilution of the aliquot with purified water and perform the assay. Multiply the result obtained by 10. If the obtained result surplus the operating range, perform a new dilution by altering the dilution proportion.

PRODUCT DESCRIPTION

R 1	Ammonium molybdate $\geq 0,5$ mmol/L; sulfuric acid $\geq 0,5\%$ v/v; triton X100 $\geq 0,1\%$ v/v.	
STD	Sodium phosphate in concentration equivalent to 5,0 mg/dL, preservative. Traceable to reference material NIST 3139a.	

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrator - Autocal H	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 5 minutes.
3. Measure the Sample's and Standard's absorbance against the blank at 340 nm. The final reaction is stable for 60 minutes.

B) CALCULATIONS

Serum:
Phosphorus (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:
Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Phosphorus (mg/dL) = Sample Absorbance x Concentration Factor

24-hour urine:
Phosphorus (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD Concentration (mg/dL)} \times 10^*$

10^* = Dilution factor applied

Phosphorus (mg/24 hours) = Phosphorus (mg/dL) * x Urinary Volume (mL)
100

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.com.br

C) INTERPRETATION

Inorganic phosphate is located in the bone as hydroxyapatite. In the extracellular fluid, phosphorus is found as inorganic phosphate and in small amounts as organic phosphate (phospholipids). In the cells, the latter is found as phospholipids, nucleic acids and adenosine triphosphate (ATP). The regulation mechanism of the levels of phosphorus and calcium ions in the blood is influenced by the interactions between parathyroid hormone and vitamin D. Phosphorus is increased in serum in cases of renal failure, hypoparathyroidism, resistance to parathyroid hormone and acromegaly. Lower levels of phosphorus occur in vitamin D deficiency, hyperparathyroidism and Fanconi's syndrome.

PERFORMANCE CHARACTERISTICS

	Operating range
Serum	0,35 a 15,00 mg/dL
Urine	6,86 a 233,85 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

	Sensitivity	
	Detection Limit	Quantification Limit
Serum	0.08 mg/dL	0.35 mg/dL
Urine	0.05 mg/dL	6.86 mg/dL

	Analytical Specificity		
	Hemoglobin	Bilirubin	Triglycerides
Serum	200 mg/dL	40 mg/dL	350 mg/dL
Urine	180 mg/dL	20 mg/dL	-

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine
Number of samples	40 in duplicate	20 in duplicate
Regression equation	$y = 0.989x + 0.069$	$y = 1.006x - 2.703$
Correlation coefficient (R)	0.9943	0.9995

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 3.0 mg/dL and 5.0 mg/dL was 1.20% and 0.28%, respectively.
Urine: by applying the regression equation, the total systematic error estimated for decision levels of 350.0 mg/24h and 900.0 mg/24h was -0.17% e 0.29%.

Precision:

Determined with two runs in duplicate per day for 20 days.

Serum Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
3.64	80	0.041	1.10	0.052	1.40
7.61	80	0.006	0.10	0.092	1.20
14.53	80	0.079	0.50	0.127	0.90

Urine Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
26,93	80	0,180	0,70	0,220	0,80
87,76	80	0,380	0,40	0,580	0,70

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	Serum	Urine
Children	4.0 – 7.0 mg/dL	400 – 1300 mg/24 hours
Adults	2.5 – 4.5 mg/dL	

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Phosphorus (mg/dL) X 0,323 = Phosphorus (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnicas products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnicas instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- YEE, H.Y. A simplified method for automated phosphorus analysis in serum and urine. *Clin. Chem.* v.14, p.898-903, 1968.
- DALY, A. J.; ERTINGSHAUSEN, G. Direct method for determining inorganic phosphate in serum with the "CentrifChem". *Clin. Chem.* v.18, p.263-265, 1972.
- GAMST, O.; TRY, K. Determination of serum phosphate without deproteinization by ultraviolet spectrophotometry of the phosphomolybdic acid complex. *Scand. J. Clin. Lab. Invest.* v.40, p.483-486, 1980.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. Tietz: Fundamentos de Química Clínica. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

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