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**ENGLISH**  
**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

**INTENDED USE**  
Kit intended to determine fructosamine in serum and plasma samples. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use. Lyophilized calibrator.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

**PREPARATION AND HANDLING**  
**Calibrator:** reconstitute with 1 mL of purified water. Homogenize gently by inversion, avoiding foam formation. Wait 30 minutes at room temperature until the product is completely dissolved.  
Reconstituted, the calibrator is stable for 15 days when stored at 2 to 8 °C, and 2 months when stored at -20 °C.

**WORKING PRINCIPLE**  
**Method:** NBT – 2<sup>nd</sup> Generation Nitrobluetetrazolium  
Fructosamine is formed by the non-enzymatic glycation of plasma proteins. In alkaline medium, fructosamine's keto-enol balance is shifted, so that the analyte can now reduce the nitrotetrazolium blue dye into a purple formazan which is spectrophotometrically measured at 550 nm. The difference between the absorbances taken within 5 and 10 minutes of incubation at 37°C is proportional to the sample's concentration of fructosamine.

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**  
**Sample Type:** serum and plasma (heparin and EDTA).  
**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.  
**Preservation:**

|                  | Temperature | Stability Period |
|------------------|-------------|------------------|
| Serum and Plasma | 4 to 8 °C   | 2 weeks          |
|                  | -20 °C      | 2 months         |

**PRODUCT DESCRIPTION**  
**R 1** Buffer ≥ 120 mmol / L, nitrotetrazolium blue (NBT) ≥ 100 µmol / L, uricase ≥ 1000 U / L, preservative, surfactants and stabilizers.

**R 2** Buffer > 1,0 mmol/L

**CAL** Lyophilized serum calibrator. Traceable to a Polylysine standard. Value printed on the bottle label\*.  
\*There are specific values for Beckman AU and other equipments. Consult the label.

**QUALITY CONTROL**  
The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the controls below:

|                                   |     |           |
|-----------------------------------|-----|-----------|
| Normal Control – Quantinorm*      | REF | 13.003.00 |
| Pathological Control – Quantialt* |     | 13.004.00 |
| Frutosamina BI                    |     | 10.036.00 |

\*Quantinorm and Quantialt controls have distinct reference ranges for Beckman AU and other equipments. Consult the control's product insert.

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 550 nm (520 – 580 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**  
**A) TEST PROCEDURE**  
1. Pipette in the test tubes:

|        | Calibrador | Sample |
|--------|------------|--------|
| CAL    | 50 µL      | -      |
| Sample | -          | 50 µL  |
| R1     | 800 µL     | 800 µL |
| R2     | 200 µL     | 200 µL |

2. Homogenize and incubate at 37 °C for 5 minutes, previously resetting the equipment to zero with purified water.
3. Measure the absorbance at 550 nm (520 – 580 nm) at 5 minutes (A<sub>1</sub>) and at 10 minutes (A<sub>2</sub>).

**B) CALCULATIONS**  
Fructosamine (µmol/L) =  $\frac{(A_2 - A_1) \text{ Sample}}{(A_2 - A_1) \text{ Calibrador}} \times \text{Calibrador value (µmol/L)}$

**Calculations with the Calibration Factor:**  
Calibration Factor =  $\frac{\text{Calibrador value (µmol/L)}}{(A_2 - A_1) \text{ Calibrador}}$

Fructosamine (µmol/L) = (A<sub>2</sub> - A<sub>1</sub>) Sample x Calibration Factor

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicalda.ind.br](http://www.biotechnicalda.ind.br)

**C) INTERPRETATION**  
Fructosamine is the generic name for ketoamines. The name refers to the structure of the ketoamine's rearrangement product formed by the interaction of glucose with the ε-amino group on the lysine residues of albumin. As glycated hemoglobin determinations, the dosage of fructosamine can be used to monitor the average blood glucose concentration over a prolonged period.

**PERFORMANCE CHARACTERISTICS,**

| Operating range        |
|------------------------|
| 25.40 to 807.00 µmol/L |

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

| Sensitivity     |                      |
|-----------------|----------------------|
| Detection Limit | Quantification Limit |
| 6.46 µmol/L     | 20.00 µmol/L         |

| Analytical Specificity |            |               |               |
|------------------------|------------|---------------|---------------|
| Bilirubin              | Hemoglobin | Triglycerides | Ascorbic Acid |
| 7.50 mg/dL             | 200 mg/dL  | 586.70 mg/dL  | 8.00 mg/dL    |
| Uric Acid              | Glucose    |               |               |
| 20 mg/dL               | 1000 mg/dL |               |               |

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

| Accuracy                    | Serum             |
|-----------------------------|-------------------|
| Number of Samples           | 40 in duplicate   |
| Regression Equation         | y = 0.987x + 1.48 |
| Correlation Coefficient (R) | 0.9941            |

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 200 µmol/L and 600 µmol/L was -0.56% and -1.05%, respectively.

**Precision:**  
Determined with two runs in duplicate per day for 20 days.

| Samples (µmol/L) | Replicates | Within-Run Precision |        | Total Precision |        |
|------------------|------------|----------------------|--------|-----------------|--------|
|                  |            | SD (µmol/L)          | CV (%) | SD (µmol/L)     | CV (%) |
| 136.43           | 80         | 1.19                 | 0.90   | 2.21            | 1.60   |
| 259.03           | 80         | 1.96                 | 0.80   | 4.04            | 1.60   |
| 478.17           | 80         | 2.60                 | 0.50   | 6.82            | 1.40   |

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

|                  |                  |
|------------------|------------------|
| Serum and Plasma | 205 - 285 µmol/L |
|------------------|------------------|

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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- ARMBRUSTER, David A. Fructosamine: Structure, Analysis, and Clinical Usefulness. *Clinical Chemistry*, v. 33, n. 12, p.2153-2163, 1987.
- BAKER, John R. et al. Use of Protein-Based Standards in Automated Colorimetric Determinations of Fructosamine in Serum. *Clinical Chemistry*, v. 31, n. 9, p.1550-1554, 1985.
- HURST, Paul L. Effect of Anticoagulants on Fructosamine Determination. *Clinical Chemistry*, v. 33, n. 10, p.1947-1948, 1987.
- KRUSE-JARRES, J. D. et al. A new colorimetric method for the determination of fructosamine. *Lab.med.*, v. 13, p.245-253, 1989.
- LIM, Yong Seng; STALEY, Michael J. Measurement of Plasma Fructosamine Evaluated for Monitoring Diabetes. *Clin. Chem.*, v. 31, nº. 5, p.731-733, 1985.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. *Tietz: Fundamentos de Química Clínica*. 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

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