## Gama GT VET

Gamma-GT VET / Gamma GT VET  
Ref. 90.016.00

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### ENGLISH

#### INTENDED USE

Kit intended to determine the activity of Gamma Glutamyl Transferase (Gamma-GT) in serum and plasma samples. Veterinary diagnostic use only.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Kinetic

GT Gamma catalyzes the transfer reaction of the glutamyl group of L-γ-glutamyl-3-carboxy-4-nitroanilide to glycylglycine, originating L-γ-glutamylglycylglycine and 5-amino-2-nitrobenzoate. The enzymatic activity is determined by the formation rate of 5-amino-2-nitrobenzoate, which can be spectrophotometrically measured at 405nm.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum and plasma EDTA.

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	30 days
	-20 °C	1 year

#### PRODUCT DESCRIPTION

**R 1** TRIS buffer ≥ 50 mmol/L; Glycylglycine ≥ 50 mmol/L; surfactant ≥ 0.5%; preservative. **X**

**R 2** L- γ -Glutamyl-3-carboxy-4-nitroaniline ≥ 10 mmol/L; MES ≥ 1 mmol/L; preservative. **X**

Traceable to the specific absorbance of the chromogen 5-amino-2-nitrobenzoate using calibrated pipettes and manual spectrophotometer that provides absolute values, based on the modified methodology from IFCC (International Federation of Clinical Chemistry and Laboratory Medicine).

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use the controls below:

Calibrator-Autocal Vet	90.039.00
Normal Control – Quantinorm Vet	90.040.00
Pathological Control – Quantialt Vet	90.041.00

**REF**

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Let the reagents reach room temperature.
2. Zero the equipment at 405 nm with purified water.
3. Pipette in the test tubes:

Sample	100 µL
R1	800 µL
Homogenize and incubate the tube for 3 minutes at 37°C	
R2	200 µL

3. Homogenize and insert it in the thermostated cuvette. Activate the chronometer.
4. After 1 minute, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

#### B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

The Gamma GT activity is calculated by multiplying the ΔA/min by the factor below:

$$\text{Gamma GT (U/L)} = \Delta A/\text{min} \times 1421$$

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

#### C) INTERPRETATION

The GT Gamma is present in the proximal renal tubule, liver, pancreas and intestine. The enzyme present in the serum originates primarily from the hepatobiliary system, however, increased GT Gamma activities occur concomitantly with the elevation of other biochemical parameters. Therefore, its analysis together with alkaline phosphatase, transaminases and bilirubins, is useful in the differential diagnosis of primary and secondary liver diseases. It is elevated in chronic alcoholics and individuals who make continuous use of anticonvulsants, such as phenytoin and phenobarbital.

### PERFORMANCE CHARACTERISTICS

Operating range
1.55 to 540.80 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
1.50 U/L

Analytical Specificity			
Species	Hemoglobin	Bilirubin	Triglycerides
Canines	200 mg/dL	18 mg/dL	750 mg/dL
Felines	200 mg/dL	30 mg/dL	500 mg/dL
Equines	500 mg/dL	18 mg/dL	1000 mg/dL
Bovines	100 mg/dL	10 mg/dL	500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Caninos	
Number of Samples	20 in duplicate
Regression Equation	y = 1.0354x – 0.2215
Correlation Coefficient (R)	0.9962
Accuracy - Felines	
Number of Samples	10 in duplicate
Regression Equation	y = 0.9796x + 0.0906
Correlation Coefficient (R)	0.9936

Accuracy - Equines	
Number of Samples	10 in duplicate
Regression Equation	y = 1.0146x + 0.1133
Correlation Coefficient (R)	0.9961
Accuracy - Bovines	
Number of Samples	10 in duplicate
Regression Equation	y = 0.9243x + 0.4102
Correlation Coefficient (R)	0.9918

#### Precisión:

Determined with two runs in duplicate.

Samples (U/L)	Replicates	Within-Run Precision		Samples (U/L)	Replicates	Total Precision	
		SD (U/L)	CV (%)			SD (U/L)	CV (%)
1,88	20	0,06	3,3	1,87	40	0,08	4,2
135,59	20	3,37	2,5	136,40	40	2,89	2,1

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Canine	1.2 – 6.4 U/L
Feline	1.3 – 5.1 U/L
Equine	4.3 – 13.4 U/L
Bovine	6.1 – 17.4 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units (µKat/L):

Gamma GT (U/L) x 0,017 = Gamma GT (µkat/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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