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ENGLISH

 **BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

INTENDED USE

Kit intended to determine the activity of Gamma Glutamyl Transferase (Gamma-GT) in serum and plasma samples. Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- The Dedicated and Semi Dedicated Reagents presentations have *onboard* stability of 30 days.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: kinetic

GT Gamma catalyzes the transfer reaction of the glutamyl group of L- γ -glutamyl-3-carboxy-4-nitroanilide to glycyglycine, originating L- γ -glutamylglycylglycine and 5-amino-2-nitrobenzoate. The enzymatic activity is determined by the formation rate of 5-amino-2-nitrobenzoate, which can be spectrophotometrically measured at 405nm.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and plasma EDTA.

Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	30 days
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1 TRIS buffer ≥ 50 mmol/L; Glycylglycine ≥ 50 mmol/L; surfactant ≥ 0.5%; preservative. 

R 2 L- γ -Glutamyl-3-carboxy-4-nitroaniline ≥ 10 mmol/L; MES ≥ 1 mmol/L; preservative. 

Traceable to the specific absorbance of the chromogen 5-amino-2-nitrobenzoate using calibrated pipettes and manual spectrophotometer that provides absolute values, based on the modified methodology from IFCC (International Federation of Clinical Chemistry and Laboratory Medicine).

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal H	REF	13.002.00
Normal Control – Quantinorm		13.003.00
Pathological Control – Quantialt		13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 405 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Let the reagents reach room temperature.
2. Zero the equipment at 405 nm with purified water.
3. Pipette in the test tubes:

R1	800 µL
Sample	100 µL
Homogenize and incubate the tube for 3 minutes at 37°C	
R2	200 µL

3. Homogenize and insert it in the thermostated cuvette. Activate the chronometer.
4. After 1 minute, measure the initial absorbance A0 and read it again every minute, for 3 minutes (A1, A2 and A3, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA /min):

$$(\Delta A/\text{min}) = \frac{(A1 - A0) + (A2 - A1) + (A3 - A2)}{3}$$

The Gamma GT activity is calculated by multiplying the ΔA /min by the factor below:

$$\text{Gamma GT (U/L)} = \Delta A/\text{min} \times 1421$$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotechnicalda.ind.br

C) INTERPRETATION

The GT Gamma is present in the proximal renal tubule, liver, pancreas and intestine. The enzyme present in the serum originates primarily from the hepatobiliary system, however, increased GT Gamma activities occur concomitantly with the elevation of other biochemical parameters. Therefore, its analysis together with alkaline phosphatase, transaminases and bilirubins, is useful in the differential diagnosis of primary and secondary liver diseases. It is elevated in chronic alcoholics and individuals who make continuous use of anticonvulsants, such as phenytoin and phenobarbital.

PERFORMANCE CHARACTERISTICS

Operating range
2,20 to 310,60 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
1,05 U/L	2,20 U/L

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
500 mg/dL	50 mg/dL	1500 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1,0571x - 0,705$
Correlation Coefficient (R)	0,9998

By applying the regression equation, the total systematic error estimated for decision levels of 40 U/L, 80 U/L and 120 U/L was 7.5%, 6.6% and 6.3%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (U/L)	Replicates	Within-Run Precision		Total Precision	
		SD (U/L)	CV (%)	SD (U/L)	CV (%)
27,16	50	1,35	5,00	1,30	4,7
113,95	50	4,33	3,18	2,62	2,3
138,54	50	6,06	4,40	4,52	3,3

CV: Coeficiente of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	37 °C
Men	< 55 U/L
Women	< 38 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (µKat/L):

Gamma GT (U/L) x 0,017 = Gamma GT (µkat/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- SCHUMANN, Gerhard et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C. Part 6. Reference Procedure for the Measurement of Catalytic Concentration of γ-Glutamyltransferase. *Clin Chem Lab Med.* v. 40 (7), p.734-738, 2002.
- SZASZ, Gabor. A Kinetic Photometric Method for Serum γ-Glutamyl Transpeptidase. *Clin. Chem.* v. 15, nº. 2, p.124-136, 1969.
- HELL, Wolfgang; EHRHARDT, Volker. Reference Ranges for Adults and Children: Pre-Analytical Considerations. *Roche Diagnostics*, p.47, 2008.
- BURTIS, Carl A.; ASHWOOD, Edward R. **Tietz: Fundamentos de Química Clínica.** 4ª. ed. Rio de Janeiro: Guanabara Koogan, 1996.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. **Tietz: Fundamentos de Química Clínica.** 6. ed. Rio de Janeiro: Saunders Elsevier, 2008.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - vol. 2**, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES		
	Consultar as instruções para utilização Consult instructions for use Consultéense las instrucciones de uso	 Descartar corretamente Dispose properly Desechar adecuadamente
REF	Número de catálogo Catalog number Número de catálogo	R <N> Reagente Reagent Reactivo
LOT	Código do lote Batch code Código de lote	 Limite de temperatura Temperature limitation Limite de temperatura
	Validade Use by date Fecha de Caducidad	X Nocivo / Irritante Harmful / Irritant Nocivo / Irritante