## Glucose Plus VET

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#### ENGLISH

#### INTENDED USE

Kit intended to determine glucose in serum and plasma samples. Veterinary diagnostic use only.


#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Enzymatic - GOD

The glucose present in the sample is oxidized by the enzyme glucose oxidase (GOD) forming gluconic acid and hydrogen peroxide. The latter oxidizes the phenolic reactant, in the presence of 4-aminoantipyrine and under the catalytic action of peroxidase, to form a pink compound (quinoneimine) that can be spectrophotometrically measured at 505 nm. The color intensity is proportional to the glucose concentration in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum and plasma (fluoride).

**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	3 days
	-20 °C	1 day

\* Samples must be immediately separated from the cellular elements.

#### PRODUCT DESCRIPTION

**R 1**

Phosphate buffer ≥ 40 mmol/L; 4-aminoantipyrine ≥ 0,2 mmol/L; phenol ≥ 1,0 mmol/L; peroxidase ≥ 500 U/L; glucose Oxidase ≥ 5000 U/L; activators; detergents and stabilizer.

**STD**

Glucose 100 mg/dL; stabilizer and preservative. Traceable to NIST reference material 917c.

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal Vet	90.039.00
Normal Control – Quantinorm Vet	90.040.00
Pathological Control – Quantialt Vet	90.041.00

**REF**

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 500 nm (450 – 550 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

- Homogenize and incubate at 37 °C for 10 minutes.
- Measure the Standard' and Sample' absorbance against the Blank at 500 nm. The final reaction is stable for 10 minutes.

##### B) CALCULATIONS

Glucose (mg/dL) =  $\frac{\text{Sample's Abs.} - \text{Blank's Abs.}}{\text{STD's Abs.} - \text{Blank's Abs.}}$  x STD Concentration (mg/dL)

##### Calculations with the Calibration Factor:

Calibration Factor =  $\frac{\text{STD Concentration (mg/dL)}}{\text{STD's Abs.} - \text{Blank's Abs.}}$

Glucose (mg/dL) = (Sample's Abs. - Blank's Abs.) x Calibration Factor

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

#### C) INTERPRETATION

Glucose is the predominant carbohydrate in the peripheral blood, being the main source of cellular energy. It is stored under the forms of glycogen in the liver or fatty acids in adipose tissue. Their amount in the blood is controlled in narrow limits by several hormones, which have the function of reducing (insulin) or increasing (glucagon, adrenaline, cortisol and hormone of the growth) its concentration. Hyperglycemia occurs in Diabetes mellitus caused by deficiency in the secretion or action of insulin. Other secondary factors can also elevate blood glucose levels, which include pancreatitis, thyroid dysfunction, renal failure or liver disease. Hypoglycemia is less prevalent, occurring in some diseases such as insulinoma, hypopituitarism or insulin-induced hypoglycemia.

#### PERFORMANCE CHARACTERISTICS

Operating range
3.83 to 741.31 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity
Quantification Limit
2.90 mg/dL

#### Analytical Specificity

Species	Hemoglobin	Bilirubin	Triglycerides
Canines	300 mg/dL	10 mg/dL	250 mg/dL
Felines	500 mg/dL	10 mg/dL	250 mg/dL
Equines	500 mg/dL	18 mg/dL	250 mg/dL
Bovines	300 mg/dL	10 mg/dL	250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	40 em duplicata
Regression Equation	y = 0,9806x + 1,1982
Correlation Coefficient (R)	0,9919

Accuracy – Felines	
Number of Samples	20 em duplicata
Regression Equation	$y = 0,9952x - 3,1139$
Correlation Coefficient (R)	0,9972
Accuracy – Equines	
Number of Samples	11 em duplicata
Regression Equation	$y = 1,0271x - 3,7932$
Correlation Coefficient (R)	0,9976
Accuracy – Bovines	
Number of Samples	10 em duplicata
Regression Equation	$y = 0,9988x - 0,4469$
Correlation Coefficient (R)	0,9947

PreCISION:

Samples (mg/dL)	Repetitions	Within-Run Precision		Samples (mg/dL)	Repetitions	Total Precision	
		SD (mg/dL)	CV (%)			SD (mg/dL)	CV (%)
70,56	20	3.19	4.5	70.99	40	2.62	3.7
289,35	20	5.85	2.0	291.75	40	8.35	2.9

Determined with two runs in duplicate.

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Canine	65 - 118 mg/dL
Feline	73 – 134 mg/dL
Equine	75 - 115 mg/dL
Bovine	45 – 75 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### Conversion to the International System of Units (mmol/L):

Glucose (mg/dL) x 0,0556 = Glucose (mmol/L)

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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