Glucose Plus

Glucose Plus / Glucosa Plus
Ref.10.046.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit intended to determine glucose in serum, plasma, urine and cerebrospinal fluid samples. Diagnostic use only.


STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Enzymatic - GOD

The glucose present in the sample is oxidized by the enzyme glucose oxidase (GOD) forming gluconic acid and hydrogen peroxide. The latter oxidizes the phenolic reactant, in the presence of 4-aminoantipyrine and under the catalytic action of peroxidase, to form a pink compound (quinoneimine) that can be spectrophotometrically measured at 505 nm. The color intensity is proportional to the glucose concentration in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum, plasma fluoride, urine and cerebrospinal fluid (CSF).
Collection and handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum, Plasma and CSF*	4 to 8 °C	3 days
	-20 °C	1 day
Urine	4 to 8 °C	2 hours
	-20 °C	2 days

* Samples must be immediately centrifuged.

PRODUCT DESCRIPTION

R 1

Phosphate buffer ≥ 40 mmol/L; 4-aminoantipyrine ≥ 0,2 mmol/L; phenol ≥ 1,0 mmol/L; peroxidase ≥ 500 U/L; glicose oxidase ≥ 5000 U/L; activators; detergents and stabilizers.

X

STD

Glucose 100 mg/dL; stabilizer and preservative. Traceable to NIST reference material 917c.

X

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the calibrator and controls below:

Calibrador – Autocal H	13.002.00
Normal Control – Quantinorm	13.003.00
Pathological Control – Quantialt	13.004.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 500 nm (450 – 550 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate at 37 °C for 10 minutes.
3. Measure the Standard's and Sample's absorbance against the Blank at 500 nm. The final reaction is stable for 10 minutes.

B) CALCULATIONS

$$\text{Glucose (mg/dL)} = \frac{\text{Sample's Abs.} - \text{Blank's Abs.}}{\text{STD's Abs.} - \text{Blank's Abs.}} \times \text{STD Concentration (mg/dL)}$$

Calculations with the Calibration Factor:

$$\text{Calibration Factor} = \frac{\text{STD Concentration (mg/dL)}}{\text{STD's Abs.} - \text{Blank's Abs.}}$$

$$\text{Glucose (mg/dL)} = (\text{Sample's Abs.} - \text{Blank's Abs.}) \times \text{Calibration Factor}$$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnica.com.br

C) INTERPRETATION

Glucose is the predominant carbohydrate in the peripheral blood, being the main source of cellular energy. its concentration in the blood is regulated by several hormones, which reduce (insulin) or increase (glucagon, adrenaline, cortisol, and hormone of the growth) its concentration. Hyperglycemia occurs mainly in Diabetes Mellitus, but can also be observed in cases of pancreatitis, thyroid dysfunction, renal and liver disease. Hypoglycemia is less prevalent, occurring in some diseases such as insulinoemia, hypopituitarism or insulin-induced hypoglycemia. In healthy individuals, glucose is not detected in urine; therefore, the elevation of glucose concentration in urine can result of the excess of blood glucose or some renal disease. In the cerebrospinal fluid, the increase in glucose concentration is associated with hyperglycemia, while its reduction is related to an increase in the Central Nervous System metabolism and to the deficient transport of glucose.

PERFORMANCE CHARACTERISTICS

Operating range	
Serum and Urina	4.5 to 800.0 mg/dL
Cerebrospinal Fluid	4.5 to 740.0 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity - Quantification Limit	
Serum	3.43 mg/dL
Urine	2.27 mg/dL
CSF	2.92 mg/dL

	Hemoglobin	Bilirubin	Triglycerides
Serum	300 mg/dL	20 mg/dL	1000 mg/dL
Urine	300 mg/dL	20 mg/dL	-
CSF	100 mg/dL	-	-

Interfering substances up to the values presented above do not cause significant alterations in the results. For CSF samples, any concentration of Transfix interferes with the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum	Urine	CSF
Number of Samples	40 in duplicates	20 in duplicates	20 in duplicates
Regression Equation	y = 1.0065x - 3.787	y = 0.9903x + 1.0367	y = 1.027x + 0.5077
Correlation Coefficient (R)	0.9998	0.9994	0.9948

Serum: by applying the obtained regression equation, the estimated systematic total error for decision levels of 45 mg/dL, 120 mg/dL and 180 mg/dL was, respectively, -7.8%, -2.5% and -1.5%.

Urine: by applying the obtained regression equation, the estimated systematic total error for decision levels of 20 mg/dL, 60 mg/dL and 100 mg/dL was, respectively, 4.2%, 0.8% and 0.1%.

CSF: by applying the obtained regression equation, the estimated systematic total error for decision levels of 40 mg/dL, 80 mg/dL and 120 mg/dL was, respectively, 4.0%, 3.3% and 3.1%.

Precision:

Serum and Urine: within-run precision was determined with 5 replicates for 10 days; between-run precision was determined in 3 conditions, with 3 replicates for 5 days. CSF: within-run precision was determined with 20 replicates for 1 day; between-run precision was determined in 2 conditions, with 20 replicates for 2 days.

Serum Samples (mg/dL)	Within-Run Precision			Serum Samples (mg/dL)			Between-Run Precision		
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)	
90.39	0.78	0.9%		89.62	0.63	0.7%			
290.82	2.16	0.7%		288.80	2.60	0.9%			
514.70	5.44	1.1%		513.45	5.40	1.1%			

Urine Samples (mg/dL)	Within-Run Precision			Urine Samples (mg/dL)			Between-Run Precision		
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)	
48.79	0.68	1.4%		45.23	0.94	1.9%			
252.91	3.12	1.2%		250.91	5.07	2.0%			
555.32	8.22	1.5%		547.59	9.66	1.8%			

CSF Samples (mg/dL)	Within-Run Precision			CSF Samples (mg/dL)			Between-Run Precision		
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)	
53.90	2.41	4.5%		53.19	2.54	4.8%			
72.72	1.53	2.1%		74.29	2.23	3.0%			
103.23	1.67	1.6%		105.14	3.45	3.3%			

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

	Serum and Plasma	Urine
Premature newborns	20 - 60 mg/dL	1 - 15 mg/dL
Children and Adults	65 - 99 mg/dL	

Cerebrospinal fluid: glucose concentrations must be 60% of the plasma concentration for samples collected at the same time.

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Glucose (mg/dL) x 0,0556 = Glucose (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotech products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotech instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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