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**ENGLISH**  
**BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.**

**INTENDED USE**  
Kit intended to determine HDL cholesterol in serum samples. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

**WORKING PRINCIPLE**  
**Method:** Phosphotungstic acid Precipitation  
The chylomicrons, the very low density lipoproteins (VLDL) and the low density lipoproteins (LDL) present in the sample are precipitated by phosphotungstic acid and magnesium ions. After centrifugation, cholesterol is determined in the supernatant containing HDL cholesterol (high density lipoproteins).

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**  
**Sample Type:** serum  
**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.  
**Preservation:**

	Temperature	Stability Period
Serum	4 to 8 °C	7 days
	-20 °C	2 months

**PRODUCT DESCRIPTION**

<b>R 1</b>	Phosphotungstic acid > 1,0 mmol / L, Magnesium chloride > 10 mmol / L.
<b>STD</b>	Butanol > 5% v / v; Triton X 100 5% v / v; stabilizer and concentration of HDL cholesterol indicated on the label.

**QUALITY CONTROL**  
The use of controls should be a routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration it is recommended to use the STD below:

HDL Cholesterol STD (included in kit)	<b>REF</b>	10.005.00
---------------------------------------	------------	-----------

- NECESSARY EQUIPMENT FOR TESTING**
- Spectrophotometer or photometer for reading at 505 nm (490 - 540 nm).
  - Thermostatic water bath at 37 °C and test tubes.
  - Glass pipettes and/or automatic, clock or chronometer.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**  
**A) TEST PROCEDURE**  
**Chylomicrons, VLDL and LDL precipitation:**  
1. Pipette in centrifuge tube:

Sample	250 µL
R1	250 µL

2. Shake vigorously for 30 seconds.  
3. Centrifuge at 3,500 rpm for 15 minutes or until a clear supernatant is obtained.  
4. Carefully collect the supernatant avoiding resuspension of the precipitate and carry out the dosage of cholesterol immediately.

**HDL Cholesterol Dosage**  
Proceed to the determination of HDL cholesterol using **Biotécnica R1 Cholesterol (CAT BT 10.004.00)** as follows:  
1. Pipette into test tubes

	Blank	STD	Sample
Purified water	100 µL	-	-
STD	-	100 µL	-
Supernatant	-	-	100 µL
Cholesterol R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize well and incubate the tubes for 10 minutes at 37 ° C.  
3. Measure the absorbance of the Standard and the Sample against the Blank at 505 nm. Color is stable for 20 minutes.

**B) CALCULATIONS**  
To correct the dilution of the samples in the precipitation (1: 2), in the calculations, multiply the value of the Standard by 2.

HDL Cholesterol (mg/dL) =  $\frac{\text{Sample Absorbance} \times [\text{Standard Concentration}(\text{mg/dL}) \times 2]}{\text{Standard Absorbance}}$

**Example:**  
Standard Concentration = 20 mg/dL  
Sample Absorbance = 0,244  
Standard Absorbance = 0,305

HDL Cholesterol (mg/dL) =  $\frac{0,244 \times [20 \times 2]}{0,305} = 32 \text{ mg/dL}$

**Calculations with the Calibration Factor:**  
Calibration Factor =  $\frac{[\text{Standard Concentration}(\text{mg/dL}) \times 2]}{\text{Standard Absorbance}}$

HDL Cholesterol (mg/dL) = Sample Absorbance x Calibration Factor

**Example:**  
Calibration Factor =  $\frac{[20 \times 2]}{0,305} = 131$

HDL Cholesterol (mg/dL) = 0,244 x 131 = 32 mg/dL

**Automation:** The HDL Cholesterol Dosing step can be applied on most automated analyzers. Instrument settings are available at [www.biotechnicaind.br](http://www.biotechnicaind.br)

**C) INTERPRETATION**  
High-density lipoproteins (HDL) are responsible for the reverse transport of cholesterol from peripheral tissues to the liver. In the liver, cholesterol is converted to bile acids which are excreted in the intestine by bile. Monitoring of HDL cholesterol is of great clinical importance, as there is an inverse correlation between serum HDL concentration and the risk of atherosclerotic disease. Reduced concentrations of HDL cholesterol, along with high triglycerides, increase cardiovascular risk.

**PERFORMANCE CHARACTERISTICS,**

<b>Operating range</b>	3.00 to 143 mg/dL
------------------------	-------------------

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
<b>Detection Limit</b>	<b>Quantification Limit</b>
0.83 mg/dL	2.03 mg/dL

Analytical Specificity		
<b>Hemoglobin</b>	<b>Bilirubin</b>	<b>Triglycerides</b>
150 mg/dL	15 mg/dL	1250 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum
Number of Samples	30 in duplicate
Regression Equation	$y = 1,031x - 0,49$
Correlation Coefficient (R)	0,9901

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 30 mg/dL and 100 mg/dL was 1,47% and 2,61%, respectively.

#### Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (µmol/L)	Replicates	Within-Run Precision		Total Precision	
		SD (µmol/L)	CV (%)	SD (µmol/L)	CV (%)
10.6	80	0.107	1.0	0.223	2.1
50.9	80	0.316	0.6	1.317	2.6
127	80	0.515	0.4	0.695	0.5

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Age Group	Desirable reference values (mg/dL)
Grown up (> 20 anos)	> 40
Children and adolescents	> 45

Reference values are desirable according to the Brazilian Consensus for the Normalization of Laboratory Determination of the Lipid Profile.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at [sac@biotechnica.ind.br](mailto:sac@biotechnica.ind.br).

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