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**ENGLISH**  
**INTENDED USE**  
Kit for the quantitative determination of human hemoglobin A1c (HbA1c). *In vitro* diagnostic use.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.
- Ready-to-use reagents. Gently homogenize before use.
- **Lyoophilized Control:** Kepp at 10-30 ° C for 10 minutes. Reconstitute with 0.5 mL of purified water. Gently homogenize by inversion avoiding foaming. After reconstitution is stable for 30 days if stored at 2 to 8 ° C

**WORKING PRINCIPLE**  
**Method:** Immunoturbidimetric  
The direct determination of HbA1c in whole blood. Total hemoglobin and HbA1c have the same non-specific absorption rate as the latex particles present in R1. When the human anti-HbA1c mouse monoclonal antibody is added (R2), the latex-HbA1c-human anti-HbA1c monoclonal antibody complex is formed. Agglutination occurs when the goat polyclonal antibody interacts with the mouse monoclonal antibody. The turbidity intensity is proportional to the concentration of HbA1c in the sample.

**SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY**  
**Sample Type:** Total Blood (EDTA).  
**Collection and handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.  
**Preservation:**

	Temperature	Stability Period
Total Blood	4 to 8 °C	7 days

**PRODUCT DESCRIPTION**

<b>R 1</b>	Latex 0,13%, buffer, stabilizer, preservative.	
<b>R 2</b>	Mouse anti-human HbA1c monoclonal antibody 0.05 mg / ml, goat anti-mouse 0.08mg / dl polyclonal IgG antibody, buffer, stabilizers, preservative.	 
<b>R HML</b>	Purified Water, preservative.	
<b>CONTROL</b>	Whole human blood, glucose and preservative. Traceable to the National Glycohemoglobin Standardization Program (NGSP). (Included in the Kit with CONTROL)	 

**QUALITY CONTROL**  
The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:  
HbA1c Calibration Set 21.009.00  
HbA1C Plus CONTROL (Included in the Kit with CONTROL) **REF** 20.025.00

**NECESSARY EQUIPMENT FOR TESTING**

- Spectrophotometer or photometer for reading at 600 - 660 nm.

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

**HEMOLYSIS**  
The procedure of hemolysis must be applied on samples, control and calibrators, as described below:

1. Hemolyse 20 µL of sample, control and calibrators previously homogenized, with 1.0 mL of RHML, in different test tubes.
2. Let the tubes rest for 5 minutes to complete hemolysis. Hemolysates are stable for up to 10 days if stored at 2 to 8 ° C.

**TEST**  
**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnica.com.br](http://www.biotechnica.com.br)

REACTION PARAMETERS	
<b>Type of reaction</b>	End Point / Ascending
<b>Wave length</b>	600 - 660 nm
<b>Temperature</b>	37 °C
<b>Sample volume/CAL*</b>	5 µL
<b>R1*</b>	180 µL
<b>Incubation time</b>	5 minutes
<b>R2*</b>	60 µL
<b>Incubation time</b>	5 minutes
<b>Calibration</b>	Calibration curve**

\* Sample volumes/CAL, R1 and R2 be modified as long as the established ratio is maintained.

**\*\*Calibration Curve:**  
Use the Calibration Set HbA1C- CAT BT 21.009.00.

Point 0	Purified water
Point 1	CAL 1
Point 2	CAL 2
Point 3	CAL 3
Point 4	CAL 4

**B) CALCULATIONS**  
The equipment automatically calculates the analyte concentration in each sample.

**C) INTERPRETATION**  
There are several subtypes of HbA1, such as HbA1a1, HbA1a2, HbA1b and HbA1c. Of these, the HbA1c fraction is the glycosylated hemoglobin itself and corresponds to 80% of the glycosylated fraction. HbA1c is formed by the binding of glucose to the valine amino acid of the N-terminal portion of the beta chain of hemoglobin by means of a stable and irreversible binding. HbA1c is found in non-diabetic adults. It reflects the average glycaemia of the last 2 to 3 months, which corresponds to the half-life of the red blood cells. The higher the glycaemia, the higher the HbA1c concentration. The diagnostic criteria for diabetes mellitus underwent several modifications, precisely due to the appearance of new evidence regarding the association of increasingly low blood glucose levels with micro and macrovascular complications, which result in retinopathy, nephropathy, neuropathy, coronary disease, cerebrovascular disease and peripheral arterial disease.

**PERFORMANCE CHARACTERISTICS**

Operating range
3 to 15%

Sensitivity
Variation of 0.087 absorbance for each 1% HbA1c, to a 4.8% level of HbA1c

Analytical Specificity	
<b>Bilirubin</b>	<b>Triglycerides</b>
50 mg/dL	2000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	34 in duplicate
Regression Equation	y = 0,996x + 0,013
Correlation Coefficient (R)	0,9995

By applying the regression equation, the total systematic error estimated for decision levels of 5.5% and 10.5% was -0.16% and -0.28%, respectively.

**Precision:**  
Determined with two runs in duplicate per day for 20 days.

Samples (%)	Replicates	Within-Run Precision		Total Precision	
		SD (%)	CV (%)	SD (%)	CV (%)
5,81	80	0,057	1,0	0,068	1,2
6,98	80	0,072	1,0	0,076	1,1
13,3	80	0,107	0,8	0,156	1,2

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs of disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

#### REFERENCE RANGES

Diagnostic criteria for diabetes mellitus recommended by the Brazilian Society of Diabetes			
Test	Normal	Pre-diabetes	Diabetes
Fasting glycemia (mg / dL)	< 100	100 a 125	≥ 126
Glucose 2 hours after oral tolerance test with 75 g of glucose (mg / dL)	< 140	140 a 199	≥ 200
HbA1c (%)	< 5,7	5,7 a 6,4	≥ 6

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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There are several subtypes of HbA1, such as HbA1a, HbA1a2, HbA1b and HbA1c. Of these, the HbA1c fraction is the glycated hemoglobin itself and corresponds to 80% of the glycated fraction. HbA1c is formed by the binding of glucose to the valine amino acid of the N-terminal portion of the beta chain of hemoglobin by means of a stable and irreversible binding. HbA1c is found in non-diabetic adults. It reflects the average glycemia of the last 2 to 3 months, which corresponds to the half-life of the red blood cells. The higher the glycemia, the higher the HbA1c concentration. The diagnostic criteria for diabetes mellitus underwent several modifications, precisely due to the appearance of new evidence regarding the association of increasingly low blood glucose levels with micro and macrovascular complications, which result in retinopathy, nephropathy, neuropathy, coronary disease, cerebrovascular disease and peripheral arterial disease.

#### CARACTERÍSTICAS DE DESEMPEÑO

Intervalo Operacional
3 a 15%

Sensibilidad
Variation of 0.087 absorbance for each 1% HbA1c, to a 4.8% level of HbA1c.

Especificidad Analítica	
Bilirrubina	Triglicéridos
50 mg/dL	2000 mg/dL

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