PERFORMANCE CHARACTERISTICS
Operating range
2.0 up to 16.0 %
Sensitivity
Quantification limit
0.40 %

Analytical Specificity		
Bilirubin	Lipids	Ascorbic Acid
≤ 50 mg/dL	≤ 700 mg/dL	≤ 50 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy was determined by comparison with a similar method by the dosage of 40 samples. The regression equation Y = 1.0572 X - 0.4156 and correlation coefficient r = 0.9909 were obtained.

Precision:

Intra-assay precision studies were performed with one batch and 20 replicates in an analytical run; inter-assay precision studies were performed with one batch and two replicates in two analytical runs.

Sample (%)	Repetitions	Within-run Precision		Total Precision	
		SD (%)	%CV	SD (%)	%CV
5.56	40	0.09	1.11	0.07	1.99
11.98	40	0.10	0.59	0.20	1.49

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.

REFERENCE RANGES

Normal	< 5.7%
Pre-diabetes	5.7 – 6.4%
Diabetes	≥ 6%

These values are for orientation only and it is recommended that each laboratory establishes its own reference range.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice, in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.

- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- Goldstein DE, Little RR, Lorenz RA, et al. Tests of glycemia in diabetes. Diabetes Care 1995;18:896-909.
- The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 1993;329:977-986.
- American Diabetes Association, "Standards of Medical Care for Patients with Diabetes Mellitus", (Position statement), Diabetes Care, 21 (Suppl.1): S23 S31 (1998).
- Pantheghini M, John WG on behalf of the IFCC Scientific Division. Implementation of haemoglobin A1c results traceable to the IFCC reference system: the way forward. Clin Chem Lab Med 2007; 45(8): 942-4.

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