+ and bind to apotransferrin. The apotransferrin-Fe3 + complex is called transferrin which carries the iron ions in the blood plasma. Approximately one third of the transferrin binding sites are occupied by Fe 3+, so the additional amount of iron that can be bound is latent iron binding capacity (CLIF). The sum of serum iron with CLIF represents the total iron binding capacity (CTIF). CTIF varies in iron metabolism disorders. In iron deficiency anemia CTIF is elevated and transferrin saturation is low. Low serum iron concentrations associated with low CTIF may occur in cases of chronic anemia, malignant tumors or infections.

PERFORMANCE CHARACTERISTICS

Operating range
39.8 to 469 µg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
17.61 µg/dL	39.8 µg/dL

Analytical Specificity		
Hemoglobin	Bilirubin	Triglycerides
150 mg/dL	20 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	Serum
Number of Samples	40 in duplicate
Regression Equation	y = 1.004x - 1.04
Correlation Coefficient (R)	0.9945

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 150 µg/dL and 300 µg/dL was -0,29% and 0.05%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (µmol/L)	Replicates	Within-Run Precision		Total Precision	
		SD (µmol/L)	CV (%)	SD (µmol/L)	CV (%)
105,5	80	1,16	1,1	3,25	3,1
173,8	80	1,16	0,7	4,45	2,6
336,2	80	2,14	0,6	4,32	1,3

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices’ instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Latent Iron Capacity - CLFF	110 to 370 µg/dL
Total Iron Capacity - CTFF	250 to 450 µg/dL
Transferrin Saturation index - IST	15 to 55 %

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): µmol/L
CTLF (µg / dL) x 0.179 = CTIF (µmol/L).

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and others current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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