Accuracy – Bovines	
Number of Samples	10 in duplicate
Regression Equation	$y = 0,9861x + 0,0298$
Correlation Coefficient (R)	0,9959

Precision:

Determined with two runs in duplicate.

Samples (mg/dL)	Within-Run Precision		Samples (mg/dL)	Total Precision	
	SD (mg/dL)	CV (%)		SD (mg/dL)	CV (%)
0,47	0,02	5,0	0,47	0,02	4,2
7,66	0,16	2,0	7,66	0,15	2,0

CV: Coefficient of variation; SD: Standard deviation

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Caninos	0,1 – 0,4 mg/dL
Felinos	0 – 0,3 mg/dL
Equinos	0,9 – 1,1 mg/dL
Bovinos	0 – 2,0 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to the International System of Units (mmol/L):

Uric Acid (mg/dL) X 59,5 = Uric Acid (mmol/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnicaltda.com.br.

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