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Rheumatic Control Level II / Control Reumático Nível II
Ref. 21.008.00

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ENGLISH

BEFORE USING THE PRODUCT, CHECK THE VERSION OF THE CORRESPONDING INSTRUCTION FOR USE ON THE LABEL.

INTENDED USE

Kit for internal quality control in assays of ASO Turbilátex (ASO), CRP Turbilátex (CRP) and RF Turbilátex (RF). Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- After reconstitution, the product is stable for 10 days when stored at 2 to 8 °C and 6 months at -20 °C.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Human lyophilized serum used to verify the performance of ASO, CRP and RF determination.

PREPARATION AND HANDLING OF THE PRODUCT

Reconstitute with 1.0 mL of purified water. Homogenize gently by inversion, avoiding foaming. Wait 10 minutes at room temperature until complete dissolution.

PRODUCT DESCRIPTION

CONTROL Level II: Human serum with concentration pathological of ASO, FR and PCR.  

The Controls present negative results for the HB antigen, for the HCV antibody. However, they should be treated as potentially infectious.

The control values are traceable to the following materials: International CRM 470 (CRP), Rheumatoid Arthritis Serum WHO 64/2

(RF) and International Standard for anti-streptolysin-O NIBSC 97/662 (ASO).

NECESSARY EQUIPMENT FOR TESTING

- Glass and/or automatic pipettes.
- ASO Turbilátex kit (BT 20.004.00).
- CRP Turbilátex Kit (BT 20.015.00).
- RF Turbilátex Kit (BT 20.007.00).
- Purified water.

TEST PROCEDURE

The controls should be used in the same way as the patient samples according to the procedure established for the ASO, CRP and RF tests and according to the equipment protocol.

PERFORMANCE CHARACTERISTICS

Precision:

It was determined using Rheumatic Control Level 2 with a daily run of 5 replicates for 5 days.

	Value	CV%
ASO	281.6 UI/mL	7.69%
FR	51.3 UI/mL	8.10%
PCR	49.5 mg/L	4.81%

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix controls from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the product if it displays any signs in disagreement with the ones specified in its MSDS.
- Avoid leaving controls outside the specified storage conditions.

WARNINGS AND PRECAUTIONS

- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at www.biotechnica.ind.br or calling +55 35 3214 4646.
- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

- All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotechnica.ind.br.
- To obtain instructions for use in printed format, at no additional cost, contact customer service: +55 35 3214 4646 or by email at sac@biotechnica.ind.br.

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- CEMBROWSKI, G. Thoughts on quality-control systems: A laboratorians perspective. Clin. Chem. v.43, p.886-892, 1997.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981. Approved Guideline. CLSI document H47-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2008

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The values reported for the control serve as a reference. Each laboratory should establish the range of variation in its own analytical system.

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