Equinos	500 mg/dL	5 mg/dL	1000 mg/dL
Bovinos	500 mg/dL	5 mg/dL	750 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in the results. For drugs, consult the recommended reference (Young, 2000).

Accuracy - Canines	
Number of Samples	20 in duplicate
Regression Equation	$y = 1.0234x - 6,5278$
Correlation Coefficient (R)	0.9975
Accuracy – Felines	
Number of Samples	13 in duplicate
Regression Equation	$y = 0.9918x + 5,5112$
Correlation Coefficient (R)	0.9986
Accuracy – Equines	
Number of Samples	10 in duplicate
Regression Equation	$y = 1.0043x + 0.0424$
Correlation Coefficient (R)	0.9999
Accuracy – Bovines	
Number of Samples	10 in duplicate
Regression Equation	$y = 0.9453x + 12.223$
Correlation Coefficient (R)	0.9975

#### Precision:

Determined with two runs in duplicate.

Samples (umol/L)	Within-Run Precision SD (umol/L)	CV (%)	Samples (umol/L)	Total Precision SD (umol/L)	CV (%)
206.05	5.2	2.5	204.75	4.48	2.2
375.24	10.02	2.7	375.45	11.06	2.9

CV: Coefficient of variation; SD: Standard deviation

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use protective equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Canines	147 - 156 mmol/L
Felines	141 - 152 mmol/L
Equines	132 - 146 mmol/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus, according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposal, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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